

Horse Resistin ELISA Kit | RESISTIN elisa kit

Horse Resistin ELISA Kit

CATALOG #
AAA10102

REACTIVITY
Horse


PRODUCT SYNONYMS

Resistin; Horse Resistin ELISA Kit; RESISTIN elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Horse
Specificity	No significant cross-reactivity or interference between Horse RESISTIN and analogues was observed.
Samples	Human Body Fluids And Tissue Homogenates
Assay Type	Quantitative Sandwich
Detection Range	1.56ng/ml-50ng/ml
Sensitivity	1.0ng/ml.
Intra-assay Precision	Intra-assay CV (%) is less than 15%
Inter-assay Precision	Inter-assay CV (%) is less than 15%. [CV(%) = SD/mean ×100]. All CV% should be compared by concentration, not compared by OD values.
Preparation and Storage	Store all reagents at 2-8 degree C

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for RESISTIN elisa kit

Background/Introduction: This Quantitative Sandwich ELISA kit is for lab reagent/research use only! This kit is intended to be used for determine the level of RSN (hereafter termed "analyte") in undiluted original Human body fluids and tissue homogenates. This kit is NOT suitable for assaying non-biological sources of substances.

NCBI AND UNIPROT INFORMATION

NCBI GI #	12584202
NCBI GeneID	56729

UniProt Accession #	Q9HD89 ; Q540D9 ; D6W649
Molecular Weight	11,419 Da
NCBI Official Full Name	resistin
NCBI Official Synonym Full Names	resistin
NCBI Official Symbol	RETN
NCBI Official Synonym Symbols	ADSF; RSTN; XCP1; FIZZ3; RETN1
NCBI Protein Information	resistin; resistin delta2; found in inflammatory zone 3; cysteine-rich secreted protein FIZZ3; adipose tissue-specific secretory factor; cysteine-rich secreted protein A12-alpha-like 2; c/EBP-epsilon-regulated myeloid-specific secreted cysteine-rich protein; C/EBP-epsilon regulated myeloid-specific secreted cysteine-rich protein precursor 1
UniProt Protein Name	Resistin
UniProt Gene Name	RETN
UniProt Synonym Gene Names	FIZZ3; HXCP1; RSTN; ADSF
UniProt Entry Name	RETN_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/horse-elisa-kits/resistin/10102>