

Mouse Vascular Endothelial cell Growth Factor D ELISA Kit | VEGFD elisa kit

Mouse Vascular Endothelial cell Growth Factor D ELISA Kit

CATALOG #
AAA101677

REACTIVITY
Mouse


PRODUCT SYNONYMS

Vascular Endothelial cell Growth Factor D; Mouse Vascular Endothelial cell Growth Factor D ELISA Kit; VEGFD/FIGF/VEGF-D/c-fos induced growth factor/vascular endothelial growth factor D/VEGF-DVEGFDc-Fos-induced growth factor; VEGFD elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Assay Type	Sandwich ELISA, Double Antibody
Samples	Serum, Plasma, Cell Culture Supernatant, cell or tissue lysate, Other liquid samples
Detection Range	31.25-2000pg/ml
Sensitivity	18.75pg/ml
Preparation and Storage	Store entire kit at 2-8C for short-term. For longer-term, please store the microplate & standard at -20C, while the remaining reagents can be stored at 2-8C

IMAGES



Typical Data & Standard Curve

Results of a typical standard run of a Mouse VEGFD ELISA Kit are shown below. This standard curve was generated at our lab for demonstration purpose only. Each user should obtain their own standard curve as per experiment (N/A-most applicable)

X	pg/ml	0	31.2	62.5	125	250	500	1000	2000
Y	OD450	0.002	0.166	0.343	0.633	1.159	1.775	1.865	2.348

Standard Curve (Sample)

Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for VEGFD elisa kit

Principle of the Assay||This kit was based on sandwich enzyme-linked immune-sorbent assay technology. Anti-Mouse VEGFD antibody was pre-coated onto 96-well plates. And the biotin conjugated anti-Mouse VEGFD antibody was used as detection antibodies. The standards, test samples and biotin conjugated detection antibody were added to the wells subsequently, and wash with wash buffer. HRPStreptavidin was added and unbound conjugates were washed away with wash buffer. TMB substrates were used to visualize HRP enzymatic reaction. TMB was catalyzed by HRP to produce a blue color product that changed into yellow after adding acidic stop solution. The density of yellow is proportional to the Mouse VEGFD amount of sample captured in plate. Read the O.D. absorbance at 450nm in a microplate reader, and then the concentration of Mouse VEGFD can be calculated.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.