

Mouse Nitric Oxide Synthase 3, Endothelial ELISA Kit | NOS3/eNOS elisa kit

Mouse Nitric Oxide Synthase 3, Endothelial ELISA Kit

CATALOG #
AAA101905

REACTIVITY
Mouse

PRODUCT SYNONYMS

Nitric Oxide Synthase 3; Endothelial; Mouse Nitric Oxide Synthase 3; Endothelial ELISA Kit; NOS3/eNOS/Constitutive NOS/ECNOS/endothelial/Endothelial NOS/eNOS/cNOS/nitric oxide synthase 3 (endothelial cell)/NOS type III/NOSIII; NOS3/eNOS elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Specificity	This assay has high sensitivity and excellent specificity for detection of NOS2/iNOS. No significant cross-reactivity or interference between NOS2/iNOS and analogues was observed.
Assay Type	Sandwich ELISA, Double Antibody
Samples	Serum, Plasma, Cell Culture Supernatant, cell or tissue lysate, Other liquid samples
Detection Range	15.625-1000pg/ml
Sensitivity	9.375pg/ml
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level NOS2/iNOS were tested 20 times on one plate, respectively. Intra-Assay: CV<8%
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level NOS2/iNOS were tested on 3 different plates, 8 replicates in each plate. CV (%) = SD/meanX100. Inter-Assay: CV<10%
Preparation and Storage	Store entire kit at 2-8C for short-term. For longer-term, please store the microplate & standard at -20C, while the remaining reagents can be stored at 2-8C

IMAGES

Typical Data & Standard Curve

Results of a typical standard operation of a NOS3/eNOS ELISA Kit are listed below. This standard curve was generated at our lab for demonstration purpose only. Users shall obtain standard curve as per experiment by themselves. (N/A-not applicable)

X	pg/ml	0	15.625	31.25	62.5	125	250	500	1000
Y	OD450	0.134	0.418	0.663	1.032	1.443	1.854	2.161	2.535

Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for NOS3/eNOS elisa kit

Principle of the Assay: This kit was based on sandwich enzyme-linked immune-sorbent assay technology. Anti-NOS2/iNOS antibody was pre-coated onto 96-well plates. And the biotin conjugated anti-NOS2/iNOS antibody was used as detection antibodies. The standards, test samples and biotin conjugated detection antibody were added to the wells subsequently, and washed with wash buffer. HRP-Streptavidin was added and unbound conjugates were washed away with wash buffer. TMB substrates were used to visualize HRP enzymatic reaction. TMB was catalyzed by HRP to produce a blue color product that changed into yellow after adding acidic stop solution. The density of yellow is proportional to the NOS2/iNOS amount of sample captured in plate. Read the O.D. absorbance at 450nm in a microplate reader, and then the concentration of NOS2/iNOS can be calculated.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

