

Human Anti-Granulocyte Macrophage Colony Stimulating Factor Antibody ELISA Kit | GM-CSF Ab elisa kit

Human Anti-Granulocyte Macrophage Colony Stimulating Factor Antibody ELISA Kit

CATALOG #
AAA102240

REACTIVITY
Human


PRODUCT SYNONYMS

Anti-Granulocyte Macrophage Colony Stimulating Factor Antibody; Human Anti-Granulocyte Macrophage Colony Stimulating Factor Antibody ELISA Kit; GM-CSF Ab; GM-CSF Ab elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of GM-CSF Ab. No significant cross-reactivity or interference between GM-CSF Ab and analogues was observed.
Assay Type	Sandwich ELISA, Double Antigen
Samples	Serum, Plasma, Cell Culture Supernatant, cell or tissue lysate, Other liquid samples
Detection Range	0.625-40ng/ml
Sensitivity	0.375ng/ml
Intra-assay Precision	Intra-Assay: CV<8%
Inter-assay Precision	Inter-Assay: CV<10%
Preparation and Storage	Store entire kit at 2-8C for short-term. For longer-term, please store the microplate & standard at -20C, while the remaining reagents can be stored at 2-8C

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for GM-CSF Ab elisa kit

Principle of the Assay: This kit was based on sandwich enzyme-linked immune-sorbent assay technology. Antigen was pre-coated onto 96-well plates. And the Biotin- labeled antigen was used as detection antigen. The standards, test samples and Biotin-labeled antigen were added to the wells subsequently, and wash with wash buffer. HRP-Streptavidin Conjugate was added and unbound conjugates were washed away with wash buffer. TMB substrates were used to visualize HRP enzymatic reaction. TMB was catalyzed by HRP to produce a blue color product that changed into yellow after adding acidic stop solution. The density of yellow is proportional to the target amount of sample captured in plate. Read the O.D. absorbance at 450nm in a microplate reader, and then the concentration of target can be calculated.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/gm-csf-ab-human-elisa-kits/anti-granulocyte-macrophage-colony-stimulating-factor-antibody/102240>