

Human Bactericidal permeability-increasing protein ELISA Kit | BPI elisa kit

Human Bactericidal permeability-increasing protein ELISA Kit

CATALOG #
AAA102392

REACTIVITY
Human


PRODUCT SYNONYMS

Bactericidal permeability-increasing protein; Human Bactericidal permeability-increasing protein ELISA Kit; BPI (Bactericidal/Permeability Increasing Protein)/Bactericidal permeability-increasing protein/CAP 57; BPI elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	487
Assay Type	Sandwich ELISA, Double Antibody
Samples	Serum, Plasma, Cell Culture Supernatant, cell or tissue lysate, Other liquid samples
Detection Range	0.313-20ng/ml
Sensitivity	0.188ng/ml
Preparation and Storage	Store entire kit at 2-8C for short-term. For longer-term, please store the microplate & standard at -20C, while the remaining reagents can be stored at 2-8C

IMAGES



Standard Curve (Sample)

Typical Data & Standard Curve

Results of a typical standard run of a Human BPI ELISA Kit are shown below. This standard curve was generated at our lab for demonstration purpose only. Each user should obtain their own standard curve as per experiment. (N/A-not applicable)

X	ng/ml	0	0.313	0.625	1.25	2.5	5	10	20
Y	OD450	0.065	0.158	0.227	0.337	0.683	1.183	1.985	2.692

Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for BPI elisa kit

Principle of the Assay||This kit was based on sandwich enzyme-linked immune-sorbent assay technology. Anti-Human BPI antibody was pre-coated onto 96-well plates. And the biotin conjugated anti-Human BPI antibody was used as detection antibodies. The standards, test samples and biotin conjugated detection antibody were added to the wells subsequently, and wash with wash buffer. HRP-Streptavidin was added and unbound conjugates were washed away with wash buffer. TMB substrates were used to visualize HRP enzymatic reaction. TMB was catalyzed by HRP to produce a blue color product that changed into yellow after adding acidic stop solution. The density of yellow is proportional to the Human BPI amount of sample captured in plate. Read the O.D. absorbance at 450nm in a microplate reader, and then the concentration of Human BPI can be calculated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	157276599
NCBI GeneID	671
NCBI Accession #	NP_001716.2
NCBI GenBank Nucleotide #	NM_001725.2
Molecular Weight	53,900 Da
NCBI Official Full Name	bactericidal permeability-increasing protein
NCBI Official Synonym Full Names	bactericidal/permeability-increasing protein
NCBI Official Symbol	BPI
NCBI Official Synonym Symbols	rBPI; BPIFD1
NCBI Protein Information	bactericidal permeability-increasing protein
UniProt Protein Name	Bactericidal permeability-increasing protein
UniProt Gene Name	BPI
UniProt Synonym Gene Names	BPI
UniProt Entry Name	BPI_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/bpi-human-elisa-kits/bactericidal-permeability-increasing-protein/102392>