

Human E Cadherin ELISA Kit | CDH1 elisa kit

Human E Cadherin ELISA Kit

CATALOG # AAA105531	REACTIVITY Human	APPLICATIONS ELISA
-------------------------------	----------------------------	------------------------------


PRODUCT SYNONYMS

E Cadherin; Human E Cadherin ELISA Kit; CDHE; UVO; Cadherin-1; CAM 120/80; Epithelial cadherin; E-cadherin; Uvomorulin; CD324; CDH1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	882
Samples	Cell Culture Supernatants, Serum, Plasma, Saliva, Urine
Assay Type	Sandwich
Detection Range	156 pg/ml - 10000 pg/ml
Sensitivity	15 pg/ml
Applicable Applications for CDH1 elisa kit	ELISA (Sandwich ELISA)
Detection Method	Colorimetric
Preparation and Storage	Shipped and store at 4 degree C for 6 months, store at -20 degree C for one year. Avoid freeze/thaw cycles.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CDH1 elisa kit	Sandwich ELISA for quantitative detection of Human E Cadherin concentrations.
--	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	4757960
NCBI GeneID	999
NCBI Accession #	NP_004351.1

NCBI GenBank Nucleotide #	NM_004360.3
UniProt Accession #	P12830; Q13799; Q14216; Q15855; Q16194; Q4PJ14; Q9UII8; A8K1U7
Molecular Weight	90,942 Da
NCBI Official Full Name	cadherin-1 preproprotein
NCBI Official Synonym Full Names	cadherin 1, type 1, E-cadherin (epithelial)
NCBI Official Symbol	CDH1
NCBI Official Synonym Symbols	UVO; CDHE; ECAD; LCAM; Arc-1; CD324
NCBI Protein Information	cadherin-1; CAM 120/80; E-Cadherin; uvomorulin; cell-CAM 120/80; epithelial cadherin; cadherin 1, E-cadherin (epithelial); calcium-dependent adhesion protein, epithelial
UniProt Protein Name	Cadherin-1
UniProt Gene Name	CDH1
UniProt Synonym Gene Names	CDHE; UVO; E-cadherin
UniProt Entry Name	CADH1_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/cdh1-human-elisa-kits/e-cadherin/105531>