

Human IL-33 ELISA Kit | IL33 elisa kit

Human IL-33 ELISA Kit

CATALOG # AAA105554	REACTIVITY Human	APPLICATIONS ELISA
-------------------------------	----------------------------	------------------------------


PRODUCT SYNONYMS

IL-33; Human IL-33 ELISA Kit; C9orf26; IL1F11; NFHEV; Interleukin-33; Interleukin-1 family member 11; IL-1F11; Nuclear factor from high endothelial venules; NF-HEV; IL33 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	228
Samples	Cell Culture Supernatants, Serum, Plasma, Cell lysates
Assay Type	Sandwich
Detection Range	15.6 pg/ml - 1000 pg/ml
Sensitivity	10 pg/ml
Applicable Applications for IL33 elisa kit	ELISA (Sandwich ELISA)
Detection Method	Colorimetric
Preparation and Storage	Shipped and store at 4 degree C for 6 months, store at -20 degree C for one year. Avoid freeze/thaw cycles.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for IL33 elisa kit	Sandwich ELISA for quantitative detection of Human IL-33 concentrations.
--	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	313851030
NCBI GeneID	90865
NCBI Accession #	NP_001186569.1

NCBI GenBank Nucleotide #	NM_001199640.1
UniProt Accession #	O95760; Q2YEJ5; B2R8L1; B4DJ35; B4E1Q9; D3DRI5; E7EAX4
Molecular Weight	16,250 Da
NCBI Official Full Name	interleukin-33 isoform 2
NCBI Official Synonym Full Names	interleukin 33
NCBI Official Symbol	IL33
NCBI Official Synonym Symbols	DVS27; IL1F11; NF-HEV; NFEHEV; C9orf26
NCBI Protein Information	interleukin-33; IL-33; IL-1F11; DVS27-related protein; interleukin-1 family member 11; nuclear factor for high endothelial venules; nuclear factor from high endothelial venules
UniProt Protein Name	Interleukin-33
UniProt Gene Name	IL33
UniProt Synonym Gene Names	C9orf26; IL1F11; NFHEV; IL-33; IL-1F11; NF-HEV
UniProt Entry Name	IL33_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/il33-human-elisa-kits/il-33/105554>