

Human CD87 ELISA Kit | PLAUR elisa kit

Human CD87 ELISA Kit

CATALOG # AAA105623	REACTIVITY Human	APPLICATIONS ELISA
-------------------------------	----------------------------	------------------------------


PRODUCT SYNONYMS

CD87; Human CD87 ELISA Kit; MO3; UPAR; Urokinase plasminogen activator surface receptor; U-PAR; uPAR; Monocyte activation antigen Mo3; PLAUR elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	281
Samples	Cell Culture Supernatants, Serum, Plasma, Urine
Assay Type	Sandwich
Detection Range	62.5 pg/ml - 4000 pg/ml
Sensitivity	4 pg/ml
Applicable Applications for PLAUR elisa kit	ELISA (Sandwich ELISA)
Detection Method	Colorimetric
Preparation and Storage	Shipped and store at 4 degree C for 6 months, store at -20 degree C for one year. Avoid freeze/thaw cycles.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PLAUR elisa kit	Sandwich ELISA for quantitative detection of Human CD87 concentrations.
---	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	53829379
NCBI GeneID	5329
NCBI Accession #	NP_001005376.1

NCBI GenBank Nucleotide #	NM_001005376.2
UniProt Accession #	Q03405; Q12876; Q15845; Q16887; Q6IB52; Q9BWT0; Q9NYC8; Q9UD69; Q9UEA6; Q9UM92; Q9UMV0; A8K409
Molecular Weight	32,016 Da
NCBI Official Full Name	urokinase plasminogen activator surface receptor isoform 2
NCBI Official Synonym Full Names	plasminogen activator, urokinase receptor
NCBI Official Symbol	PLAUR
NCBI Official Synonym Symbols	CD87; UPAR; URKR; U-PAR
NCBI Protein Information	urokinase plasminogen activator surface receptor; monocyte activation antigen Mo3; u-plasminogen activator receptor form 2; urokinase-type plasminogen activator (uPA) receptor
UniProt Protein Name	Urokinase plasminogen activator surface receptor
UniProt Gene Name	PLAUR
UniProt Synonym Gene Names	MO3; UPAR; U-PAR; uPAR
UniProt Entry Name	UPAR_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/plaur-human-elisa-kits/cd87/105623>