

Mouse VEGFR2 ELISA Kit | KDR elisa kit

Mouse VEGFR2 ELISA Kit

CATALOG # AAA105641	REACTIVITY Mouse	APPLICATIONS ELISA
-------------------------------	----------------------------	------------------------------


PRODUCT SYNONYMS

VEGFR2; Mouse VEGFR2 ELISA Kit; FLK1; Vascular endothelial growth factor receptor 2; VEGFR-2; Fetal liver kinase 1; FLK-1; Kinase insert domain receptor; KDR; Protein-tyrosine kinase receptor flk-1; CD309; KDR elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Sequence Length	1367
Samples	Cell Culture Supernatants, Serum
Assay Type	Sandwich
Detection Range	62.5 pg/ml - 4000 pg/ml
Sensitivity	2 pg/ml
Applicable Applications for KDR elisa kit	ELISA (Sandwich ELISA)
Detection Method	Colorimetric
Preparation and Storage	Shipped and store at 4 degree C for 6 months, store at -20 degree C for one year. Avoid freeze/thaw cycles.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for KDR elisa kit	Sandwich ELISA for quantitative detection of Mouse VEGFR2 concentrations.
---	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	549321
NCBI GeneID	16542
UniProt Accession #	P35918; C5H7S5

Molecular Weight	75,230 Da
NCBI Official Full Name	Vascular endothelial growth factor receptor 2
NCBI Official Synonym Full Names	kinase insert domain protein receptor
NCBI Official Symbol	Kdr
NCBI Official Synonym Symbols	Flk1; Ly73; Flk-1; Krd-1; VEGFR2; VEGFR-2; sVEGFR-2; 6130401C07
NCBI Protein Information	vascular endothelial growth factor receptor 2; kinase NYK; VEGF receptor-2; fetal liver kinase 1; protein-tyrosine kinase receptor flk-1; vascular endothelial growth factor receptor-2; vascular endothelial growth factor receptor-3; vascular endothelial growth factor receptor- 2; soluble vascular endothelial growth factor receptor 2
UniProt Protein Name	Vascular endothelial growth factor receptor 2
UniProt Gene Name	Kdr
UniProt Synonym Gene Names	Flk-1; Flk1; VEGFR-2; FLK-1
UniProt Entry Name	VGFR2_MOUSE

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/kdr-mouse-elisa-kits/vegfr2/105641>