

Human CD244 ELISA Kit | CD244 elisa kit

Human CD244 ELISA Kit

CATALOG # AAA105688	REACTIVITY Human	APPLICATIONS ELISA
-------------------------------	----------------------------	------------------------------


PRODUCT SYNONYMS

CD244; Human CD244 ELISA Kit; 2B4; Natural killer cell receptor 2B4; NK cell activation-inducing ligand; NAIL; NK cell type I receptor protein 2B4; NKR2B4; h2B4; SLAM family member 4; SLAMF4; Signaling lymphocytic activation molecule 4; CD244 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	370
Samples	Cell Culture Supernatants, Serum, Plasma
Assay Type	Sandwich
Detection Range	62.5 pg/ml - 4000 pg/ml
Sensitivity	10 pg/ml
Applicable Applications for CD244 elisa kit	ELISA (Sandwich ELISA)
Detection Method	Colorimetric
Preparation and Storage	Shipped and store at 4 degree C for 6 months, store at -20 degree C for one year. Avoid freeze/thaw cycles.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CD244 elisa kit	Sandwich ELISA for quantitative detection of Human CD244 concentrations.
---	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	262263435
NCBI GeneID	51744
NCBI Accession #	NP_001160135.1

NCBI GenBank Nucleotide #	NM_001166663.1
UniProt Accession #	Q9BZW8; Q5VYI2; Q5VYI6; Q5VYI7; Q96T47; Q9NQD2; Q9NQD3; Q9Y288
Molecular Weight	30,784 Da
NCBI Official Full Name	natural killer cell receptor 2B4 isoform 2
NCBI Official Synonym Full Names	CD244 molecule, natural killer cell receptor 2B4
NCBI Official Symbol	CD244
NCBI Official Synonym Symbols	2B4; NAIL; Nmrk; NKR2B4; SLAMF4
NCBI Protein Information	natural killer cell receptor 2B4; h2B4; SLAM family member 4; NK cell activation-inducing ligand; NK cell type I receptor protein 2B4; NK cell activation inducing ligand NAIL; signaling lymphocytic activation molecule 4
UniProt Protein Name	Natural killer cell receptor 2B4
UniProt Gene Name	CD244
UniProt Synonym Gene Names	2B4; NAIL; NKR2B4; h2B4; SLAMF4
UniProt Entry Name	CD244_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/cd244/105688>