

Human Uteroglobulin ELISA Kit | SCGB1A1 elisa kit

Human Uteroglobulin ELISA Kit

CATALOG # AAA105766	REACTIVITY Human	APPLICATIONS ELISA
-------------------------------	----------------------------	------------------------------


PRODUCT SYNONYMS

Uteroglobulin; Human Uteroglobulin ELISA Kit; CC10; CCSP; UGB; Clara cell phospholipid-binding protein; CCPBP; Clara cells 10 kDa secretory protein; Secretoglobulin family 1A member 1; Urinary protein 1; UP-1; UP1; Urine protein 1; SCGB1A1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	91
Samples	Cell Culture Supernatants, Serum, Plasma, Saliva, Urine
Assay Type	Sandwich
Detection Range	156 pg/ml - 10000 pg/ml
Sensitivity	10 pg/ml
Applicable Applications for SCGB1A1 elisa kit	ELISA (Sandwich ELISA)
Detection Method	Colorimetric
Preparation and Storage	Shipped and store at 4 degree C for 6 months, store at -20 degree C for one year. Avoid freeze/thaw cycles.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for SCGB1A1 elisa kit	Sandwich ELISA for quantitative detection of Human Uteroglobulin concentrations.
---	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	4507809
NCBI GeneID	7356
NCBI Accession #	NP_003348.1

NCBI GenBank Nucleotide #	NM_003357.4
UniProt Accession #	P11684 ; Q6FHH3 ; Q9UCM2 ; Q9UCM4 ; B2R5F2
Molecular Weight	9,994 Da
NCBI Official Full Name	uteroglobin
NCBI Official Synonym Full Names	secretoglobin, family 1A, member 1 (uteroglobin)
NCBI Official Symbol	SCGB1A1
NCBI Official Synonym Symbols	UGB; UP1; CC10; CC16; CCSP; UP-1; CCPBP
NCBI Protein Information	uteroglobin; blastokinin; urine protein 1; urinary protein 1; clara cells 10 kDa secretory protein; clara cell phospholipid-binding protein
UniProt Protein Name	Uteroglobin
UniProt Gene Name	SCGB1A1
UniProt Synonym Gene Names	CC10; CCSP; UGB; CCPBP; CC10; UP-1; UP1; Urine protein 1
UniProt Entry Name	UTER_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/scgb1a1-human-elisa-kits/uteroglobin/105766>