

Human Ficolin 3 ELISA Kit | FCN3 elisa kit

Human Ficolin 3 ELISA Kit

CATALOG # AAA105773	REACTIVITY Human	APPLICATIONS ELISA
-------------------------------	----------------------------	------------------------------


PRODUCT SYNONYMS

Ficolin 3; Human Ficolin 3 ELISA Kit; FCNH; HAKA1; Ficolin-3; Collagen/fibrinogen domain-containing lectin 3 p35; Collagen/fibrinogen domain-containing protein 3; Hakata antigen; FCN3 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	299
Samples	Cell Culture Supernatants, Serum, Plasma
Assay Type	Sandwich
Detection Range	780 pg/ml - 50000 pg/ml
Sensitivity	20 pg/ml
Applicable Applications for FCN3 elisa kit	ELISA (Sandwich ELISA)
Detection Method	Colorimetric
Preparation and Storage	Shipped and store at 4 degree C for 6 months, store at -20 degree C for one year. Avoid freeze/thaw cycles.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for FCN3 elisa kit	Sandwich ELISA for quantitative detection of Human Ficolin 3 concentrations.
--	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	27754776
NCBI GeneID	8547
NCBI Accession #	NP_003656.2

NCBI GenBank Nucleotide #	NM_003665.2
UniProt Accession #	O75636; Q6IBJ5; Q8WW86
Molecular Weight	31,678 Da
NCBI Official Full Name	ficolin-3 isoform 1
NCBI Official Synonym Full Names	ficolin (collagen/fibrinogen domain containing) 3
NCBI Official Symbol	FCN3
NCBI Official Synonym Symbols	FCNH; HAKA1
NCBI Protein Information	ficolin-3; H-ficolin; ficolin 3; collagen/fibrinogen domain-containing protein 3; collagen/fibrinogen domain-containing lectin 3 p35; ficolin (collagen/fibrinogen domain containing) 3 (Hakata antigen)
UniProt Protein Name	Ficolin-3
UniProt Gene Name	FCN3
UniProt Synonym Gene Names	FCNH; HAKA1
UniProt Entry Name	FCN3_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/fcn3-human-elisa-kits/ficolin-3/105773>