

Human CD62P ELISA Kit | SELP elisa kit

Human CD62P ELISA Kit

CATALOG # AAA105799	REACTIVITY Human	APPLICATIONS ELISA
-------------------------------	----------------------------	------------------------------


PRODUCT SYNONYMS

CD62P; Human CD62P ELISA Kit; GMRP; GRMP; P-selectin; CD62 antigen-like family member P; Granule membrane protein 140; GMP-140; Leukocyte-endothelial cell adhesion molecule 3; LECAM3; Platelet activation dependent granule-external membrane protein; PADGEM; SELP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	830
Samples	Cell Culture Supernatants, Serum, Plasma
Assay Type	Sandwich
Detection Range	62.5 pg/ml - 4000 pg/ml
Sensitivity	5 pg/ml
Applicable Applications for SELP elisa kit	ELISA (Sandwich ELISA)
Detection Method	Colorimetric
Preparation and Storage	Shipped and store at 4 degree C for 6 months, store at -20 degree C for one year. Avoid freeze/thaw cycles.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for SELP elisa kit	Sandwich ELISA for quantitative detection of Human CD62P concentrations.
--	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	157419154
NCBI GeneID	6403
NCBI Accession #	NP_002996.2

NCBI GenBank Nucleotide #	NM_003005.3
UniProt Accession #	P16109; Q5R344; Q8IVD1
Molecular Weight	90,834 Da
NCBI Official Full Name	P-selectin
NCBI Official Synonym Full Names	selectin P (granule membrane protein 140kDa, antigen CD62)
NCBI Official Symbol	SELP
NCBI Official Synonym Symbols	CD62; GRMP; PSEL; CD62P; GMP140; LECAM3; PADGEM
NCBI Protein Information	P-selectin; GMP-140; granule membrane protein 140; granulocyte membrane protein; CD62 antigen-like family member P; platelet alpha-granule membrane protein; leukocyte-endothelial cell adhesion molecule 3; platelet activation dependent granule-external membrane protein
UniProt Protein Name	P-selectin
UniProt Gene Name	SELP
UniProt Synonym Gene Names	GMRP; GRMP; GMP-140; LECAM3; PADGEM
UniProt Entry Name	LYAM3_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/selp-human-elisa-kits/cd62p/105799>