

Human CD143 ELISA Kit | ACE elisa kit

Human CD143 ELISA Kit

CATALOG # AAA105852	REACTIVITY Human	APPLICATIONS ELISA
-------------------------------	----------------------------	------------------------------


PRODUCT SYNONYMS

CD143; Human CD143 ELISA Kit; DCP; DCP1; Angiotensin-converting enzyme; ACE; Dipeptidyl carboxypeptidase I; Kininase II; ACE elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	1306
Samples	Cell Culture Supernatants, Serum, Plasma, Saliva
Assay Type	Sandwich
Detection Range	156 pg/ml - 10000 pg/ml
Sensitivity	5 pg/ml
Applicable Applications for ACE elisa kit	ELISA (Sandwich ELISA)
Detection Method	Colorimetric
Preparation and Storage	Shipped and store at 4 degree C for 6 months, store at -20 degree C for one year. Avoid freeze/thaw cycles.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for ACE elisa kit	Sandwich ELISA for quantitative detection of Human CD143 concentrations.
---	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	4503273
NCBI GeneID	1636
NCBI Accession #	NP_000780.1

NCBI GenBank Nucleotide #	NM_000789.3
UniProt Accession #	P12821; P22966; Q53YX9; Q59GY8; Q7M4L4; B0LPF0; B4DXI3; E7EU16
Molecular Weight	78,694 Da
NCBI Official Full Name	angiotensin-converting enzyme isoform 1
NCBI Official Synonym Full Names	angiotensin I converting enzyme
NCBI Official Symbol	ACE
NCBI Official Synonym Symbols	DCP; ICH; ACE1; DCP1; CD143; MVCD3
NCBI Protein Information	angiotensin-converting enzyme; kininase II; peptidase P; CD143 antigen; testicular ECA; carboxycathepsin; dipeptidyl carboxypeptidase 1; dipeptidyl carboxypeptidase I; angiotensin converting enzyme, somatic isoform; angiotensin I converting enzyme (peptidyl-dipeptidase A) 1; angiotensin I converting enzyme peptidyl-dipeptidase A 1 transcript
UniProt Protein Name	Angiotensin-converting enzyme
UniProt Gene Name	ACE
UniProt Synonym Gene Names	DCP; DCP1; ACE
UniProt Entry Name	ACE_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/ace-human-elisa-kits/cd143/105852>