

Human FAP alpha ELISA Kit | FAP elisa kit

Human FAP alpha ELISA Kit

CATALOG # AAA105964	REACTIVITY Human	APPLICATIONS ELISA
-------------------------------	----------------------------	------------------------------


PRODUCT SYNONYMS

FAP alpha; Human FAP alpha ELISA Kit; Seprase; 170 kDa melanoma membrane-bound gelatinase; Fibroblast activation protein alpha; Integral membrane serine protease; FAP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	735
Samples	Cell Culture Supernatants, Serum, Plasma
Assay Type	Sandwich
Detection Range	62.5 pg/ml - 4000 pg/ml
Sensitivity	10 pg/ml
Applicable Applications for FAP elisa kit	ELISA (Sandwich ELISA)
Detection Method	Colorimetric
Preparation and Storage	Shipped and store at 4 degree C for 6 months, store at -20 degree C for one year. Avoid freeze/thaw cycles.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for FAP elisa kit	Sandwich ELISA for quantitative detection of Human FAP alpha concentrations.
---	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	620597362
NCBI GeneID	2191
NCBI Accession #	NP_001278736.1

NCBI GenBank Nucleotide #	NM_001291807.1
UniProt Accession #	Q12884 ; O00199 ; Q53TP5 ; Q86Z29 ; Q99998 ; Q9UID4
Molecular Weight	26,954 Da
NCBI Official Full Name	prolyl endopeptidase FAP isoform 2
NCBI Official Synonym Full Names	fibroblast activation protein, alpha
NCBI Official Symbol	FAP
NCBI Official Synonym Symbols	FAPA; SIMP; DPPIV
NCBI Protein Information	prolyl endopeptidase FAP; seprase; FAPalpha; dipeptidyl peptidase FAP; surface-expressed protease; post-proline cleaving enzyme; gelatine degradation protease FAP; integral membrane serine protease; serine integral membrane protease; 170 kDa melanoma membrane-bound gelatinase
UniProt Protein Name	Prolyl endopeptidase FAP
UniProt Gene Name	FAP
UniProt Entry Name	SEPR_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/fap-alpha/105964>