

Human Serpin A4 ELISA Kit | SERPINA4 elisa kit

Human Serpin A4 ELISA Kit

CATALOG # AAA105982	REACTIVITY Human	APPLICATIONS ELISA
-------------------------------	----------------------------	------------------------------


PRODUCT SYNONYMS

Serpine A4; Human Serpin A4 ELISA Kit; KST; PI4; Kallistatin; Kallikrein inhibitor; Peptidase inhibitor 4; PI-4; SERPINA4 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	464
Samples	Cell Culture Supernatants, Serum, Plasma
Assay Type	Sandwich
Detection Range	156 pg/ml - 10000 pg/ml
Sensitivity	10 pg/ml
Applicable Applications for SERPINA4 elisa kit	ELISA (Sandwich ELISA)
Detection Method	Colorimetric
Preparation and Storage	Shipped and store at 4 degree C for 6 months, store at -20 degree C for one year. Avoid freeze/thaw cycles.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for SERPINA4 elisa kit	Sandwich ELISA for quantitative detection of Human Serpin A4 concentrations.
--	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	573459706
NCBI GeneID	5267
NCBI Accession #	NP_001275961.1

NCBI GenBank Nucleotide #	NM_001289032.1
UniProt Accession #	P29622 ; Q53XB5 ; Q86TR9 ; Q96BZ5
Molecular Weight	48,542 Da
NCBI Official Full Name	kallistatin isoform 1
NCBI Official Synonym Full Names	serpin peptidase inhibitor, clade A (alpha-1 antiproteinase, antitrypsin), member 4
NCBI Official Symbol	SERPINA4
NCBI Official Synonym Symbols	KAL; KST; PI4; KLST; PI-4; kallistatin
NCBI Protein Information	kallistatin; serpin A4; kallikrein inhibitor; peptidase inhibitor 4; protease inhibitor 4 (kallistatin); serine (or cysteine) proteinase inhibitor, clade A (alpha-1 antiproteinase, antitrypsin), member 4
UniProt Protein Name	Kallistatin
UniProt Gene Name	SERPINA4
UniProt Synonym Gene Names	KST; PI4; PI-4
UniProt Entry Name	KAIN_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/serpina4-human-elisa-kits/serpin-a4/105982>