

Human CD40 ELISA Kit | CD40 elisa kit

Human CD40 ELISA Kit

CATALOG # AAA106002	REACTIVITY Human	APPLICATIONS ELISA
-------------------------------	----------------------------	------------------------------


PRODUCT SYNONYMS

CD40; Human CD40 ELISA Kit; TNFRSF5; Tumor necrosis factor receptor superfamily member 5; B-cell surface antigen CD40; Bp50; CD40L receptor; CDw40; CD40 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	277
Samples	Cell Culture Supernatants, Serum, Plasma, Saliva, Urine
Assay Type	Sandwich
Detection Range	15.6 pg/ml - 1000 pg/ml
Sensitivity	1 pg/ml
Applicable Applications for CD40 elisa kit	ELISA (Sandwich ELISA)
Detection Method	Colorimetric
Preparation and Storage	Shipped and store at 4 degree C for 6 months, store at -20 degree C for one year. Avoid freeze/thaw cycles.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CD40 elisa kit	Sandwich ELISA for quantitative detection of Human CD40 concentrations.
--	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	4507581
NCBI GeneID	958
NCBI Accession #	NP_001241.1

NCBI GenBank Nucleotide #	NM_001250.5
UniProt Accession #	P25942; Q53GN5; Q5JY15; Q5U007; Q7M4Q8; Q86YK5; Q9BYU0; E1P5S9
Molecular Weight	22,259 Da
NCBI Official Full Name	tumor necrosis factor receptor superfamily member 5 isoform 1
NCBI Official Synonym Full Names	CD40 molecule, TNF receptor superfamily member 5
NCBI Official Symbol	CD40
NCBI Official Synonym Symbols	p50; Bp50; CDW40; TNFRSF5
NCBI Protein Information	tumor necrosis factor receptor superfamily member 5; CD40L receptor; CD40 type II isoform; B cell-associated molecule; B cell surface antigen CD40; CD40 antigen (TNF receptor superfamily member 5); nerve growth factor receptor-related B-lymphocyte activation molecule
UniProt Protein Name	Tumor necrosis factor receptor superfamily member 5
UniProt Gene Name	CD40
UniProt Synonym Gene Names	TNFRSF5
UniProt Entry Name	TNR5_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
 Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/cd40/106002>