

Human Ficolin 2 ELISA Kit | FCN2 elisa kit

Human Ficolin 2 ELISA Kit

CATALOG # AAA106008	REACTIVITY Human	APPLICATIONS ELISA
-------------------------------	----------------------------	------------------------------


PRODUCT SYNONYMS

Ficolin 2; Human Ficolin 2 ELISA Kit; FCNL; Ficolin-2; 37 kDa elastin-binding protein; Collagen/fibrinogen domain-containing protein 2; EBP-37; Ficolin-B; Ficolin-beta; Hucolin; L-ficolin; Serum lectin p35; FCN2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	313
Samples	Cell Culture Supernatants, Serum, Plasma
Assay Type	Sandwich
Detection Range	780 pg/ml - 50000 pg/ml
Sensitivity	20 pg/ml
Applicable Applications for FCN2 elisa kit	ELISA (Sandwich ELISA)
Detection Method	Colorimetric
Preparation and Storage	Shipped and store at 4 degree C for 6 months, store at -20 degree C for one year. Avoid freeze/thaw cycles.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for FCN2 elisa kit	Sandwich ELISA for quantitative detection of Human Ficolin 2 concentrations.
--	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	61744445
NCBI GeneID	2220
NCBI Accession #	NP_004099.2

NCBI GenBank Nucleotide #	NM_004108.2
UniProt Accession #	Q15485; Q6IS69; Q7M4P4; Q9UC57; A6NFG7; A8K478
Molecular Weight	30,228 Da
NCBI Official Full Name	ficolin-2 isoform a
NCBI Official Synonym Full Names	ficolin (collagen/fibrinogen domain containing lectin) 2
NCBI Official Symbol	FCN2
NCBI Official Synonym Symbols	P35; FCNL; EBP-37; ficolin-2
NCBI Protein Information	ficolin-2; L-ficolin; ficolin B; ficolin-B; ficolin-beta; serum lectin p35; 37 kDa elastin-binding protein; collagen/fibrinogen domain-containing protein 2; ficolin (collagen/fibrinogen domain containing lectin) 2 (hucolin)
UniProt Protein Name	Ficolin-2
UniProt Gene Name	FCN2
UniProt Synonym Gene Names	FCNL
UniProt Entry Name	FCN2_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/fcn2-human-elisa-kits/ficolin-2/106008>