

Rat PDGF-AB ELISA Kit | PDGFA elisa kit

Rat PDGF-AB ELISA Kit

CATALOG # AAA106070	REACTIVITY Rat	APPLICATIONS ELISA
-------------------------------	--------------------------	------------------------------

PRODUCT SYNONYMS

PDGF-AB; Rat PDGF-AB ELISA Kit; PDGF1; Platelet-derived growth factor subunit A; PDGF subunit A; PDGF-1; Platelet-derived growth factor A chain; Platelet-derived growth factor alpha polypeptide; PDGF2; SIS; Platelet-derived growth factor subunit B; PDGF subunit B; PDGF-2; Platelet-derived growth factor B chain; Platelet-derived growth factor beta polypeptide; Proto-oncogene c-Sis; Becaplermin; PDGFA elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Rat
Sequence Length	196
Samples	Cell Culture Supernatants, Serum, Plasma
Assay Type	Sandwich
Detection Range	31.2 pg/ml - 2000 pg/ml
Sensitivity	3 pg/ml
Applicable Applications for PDGFA elisa kit	ELISA (Sandwich ELISA)
Detection Method	Colorimetric
Preparation and Storage	Shipped and store at 4 degree C for 6 months, store at -20 degree C for one year. Avoid freeze/thaw cycles.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PDGFA elisa kit	Sandwich ELISA for quantitative detection of Rat PDGF-AB concentrations.
---	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	29789040
NCBI GeneID	25266
NCBI Accession #	NP_036933.1
NCBI GenBank Nucleotide #	NM_012801.1
UniProt Accession #	P28576; Q05028

Molecular Weight	22,338 Da
NCBI Official Full Name	platelet-derived growth factor subunit A
NCBI Official Synonym Full Names	platelet-derived growth factor alpha polypeptide
NCBI Official Symbol	Pdgfa
NCBI Official Synonym Symbols	PDGFACP
NCBI Protein Information	platelet-derived growth factor subunit A; PDGF-1; PDGF subunit A; platelet derived growth factor, alpha; Platelet-derived growth factor A chain
UniProt Protein Name	Platelet-derived growth factor subunit A
UniProt Gene Name	Pdgfa
UniProt Synonym Gene Names	Rpa1; PDGF subunit A
UniProt Entry Name	PDGFA_RAT

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/pdgfa-rat-elisa-kits/pdgf-ab/106070>