

Human calpain 5 ELISA Kit | CAPN5 elisa kit

Human Calpain-5, CAPN5 ELISA Kit

CATALOG #
AAA112549

REACTIVITY
Human


PRODUCT SYNONYMS

calpain 5; Human Calpain-5; CAPN5 ELISA Kit; Human Calpain-5 (CAPN5) ELISA kit; FLJ46245; HTRA3; nCL-3; ; CAPN5 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of human CAPN5. No significant cross-reactivity or interference between human CAPN5 and analogues was observed.
Sequence Length	640
Samples	Serum, plasma, tissue homogenates, Cell lysates
Assay Type	Sandwich
Detection Range	31.25 pg/ml -2000 pg/ml.
Sensitivity	The minimum detectable dose of human CAPN5 is typically less than 7.8 pg/ml. The sensitivity of this assay, or Lower Limit of Detection (LLD) was defined as the lowest protein concentration that could be differentiated from zero. It was determined the mean O.D value of 20 replicates of the zero standard added by their three standard deviations.
Intra-assay Precision	Intra-assay Precision (Precision within an assay): CV%<8%. Three samples of known concentration were tested twenty times on one plate to assess.
Inter-assay Precision	Inter-assay Precision (Precision between assays): CV%<10%. Three samples of known concentration were tested in twenty assays to assess.
Preparation and Storage	Unopened test kits should be stored at 2 to 8 degree C upon receipt. Please refer to pdf manual for further storage instructions.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CAPN5 elisa kit	Principle of the assay: This assay employs the quantitative sandwich enzyme immunoassay technique. Antibody specific for CAPN5 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any CAPN5 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for CAPN5 is added to the wells. After washing, avidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of CAPN5 bound in the initial step. The color development is stopped and the intensity of the color is measured.
---	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	37577157
NCBI GeneID	726
NCBI Accession #	NP_004046.2
NCBI GenBank Nucleotide #	NM_004055.4
UniProt Accession #	O15484; O00263
Molecular Weight	73,169 Da
NCBI Official Full Name	calpain-5
NCBI Official Synonym Full Names	calpain 5
NCBI Official Symbol	CAPN5
NCBI Official Synonym Symbols	VRNI; ADNIV; HTRA3; nCL-3
NCBI Protein Information	calpain-5; new calpain 3; calpain htra-3
UniProt Protein Name	Calpain-5
UniProt Gene Name	CAPN5
UniProt Synonym Gene Names	NCL3; nCL-3
UniProt Entry Name	CAN5_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/capn5-human-elisa-kits/calpain-5/112549>