

Human peroxiredoxin 1 ELISA Kit | PRDX1 elisa kit

Human Peroxiredoxin-1, PRDX1 ELISA Kit

CATALOG #
AAA112740

REACTIVITY
Human


PRODUCT SYNONYMS

peroxiredoxin 1; Human Peroxiredoxin-1; PRDX1 ELISA Kit; Human Peroxiredoxin-1 (PRDX1) ELISA kit; MSP23; NKEFA; PAG; PAGA; PAGB; PRX1; PRXI; TDPX2; natural killer-enhancing factor A; proliferation-associated gene A; thioredoxin-dependent peroxide reductase 2; PRDX1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of human PRDX1. No significant cross-reactivity or interference between human PRDX1 and analogues was observed.
Sequence Length	199
Samples	Serum, plasma, tissue homogenates, Cell lysates
Assay Type	Sandwich
Detection Range	3.12 ng/ml -200 ng/ml.
Sensitivity	The minimum detectable dose of human PRDX1 is typically less than 0.78 ng/ml. The sensitivity of this assay, or Lower Limit of Detection (LLD) was defined as the lowest protein concentration that could be differentiated from zero. It was determined the mean O.D value of 20 replicates of the zero standard added by their three standard deviations.
Intra-assay Precision	Intra-assay Precision (Precision within an assay): CV%<8%. Three samples of known concentration were tested twenty times on one plate to assess.
Inter-assay Precision	Inter-assay Precision (Precision between assays): CV%<10%. Three samples of known concentration were tested in twenty assays to assess.
Preparation and Storage	Unopened test kits should be stored at 2 to 8 degree C upon receipt. Please refer to pdf manual for further storage instructions.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PRDX1 elisa kit	Principle of the assay: This assay employs the quantitative sandwich enzyme immunoassay technique. Antibody specific for PRDX1 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any PRDX1 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for PRDX1 is added to the wells. After washing, avidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of PRDX1 bound in the initial step. The color development is stopped and the intensity of the color is measured.
---	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	320461711
NCBI GeneID	5052
NCBI Accession #	NP_001189360.1
NCBI GenBank Nucleotide #	NM_001202431.1
UniProt Accession #	Q06830; P35703; Q2V576; Q5T154; Q5T155; B5BU26; D3DPZ8
Molecular Weight	22,110 Da
NCBI Official Full Name	peroxiredoxin-1
NCBI Official Synonym Full Names	peroxiredoxin 1
NCBI Official Symbol	PRDX1
NCBI Official Synonym Symbols	PAG; PAGA; PAGB; PRX1; PRXI; MSP23; NKEFA; TDPX2; NKEF-A
NCBI Protein Information	peroxiredoxin-1; thioredoxin peroxidase 2; proliferation-associated gene A; natural killer-enhancing factor A; proliferation-associated gene protein; natural killer cell-enhancing factor A; thioredoxin-dependent peroxide reductase 2
UniProt Protein Name	Peroxiredoxin-1
UniProt Gene Name	PRDX1
UniProt Synonym Gene Names	PAGA; PAGB; TDPX2; NKEF-A; PAG
UniProt Entry Name	PRDX1_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/prdx1-human-elisa-kits/peroxiredoxin-1/112740>