

Human kallikrein-related peptidase 5 ELISA Kit | KLK5 elisa kit

Human Kallikrein-5, KLK5 ELISA Kit

CATALOG #
AAA112782

REACTIVITY
Human


PRODUCT SYNONYMS

kallikrein-related peptidase 5; Human Kallikrein-5; KLK5 ELISA Kit; Human Kallikrein-5 (KLK5) ELISA kit; KLK-L2; KLKL2; SCTE; kallikrein 5; kallikrein-like protein 2; stratum corneum tryptic enzyme; KLK5 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of human KLK5. No significant cross-reactivity or interference between human KLK5 and analogues was observed.
Sequence Length	293
Samples	Serum, plasma, cell culture supernates, tissue homogenates, saliva
Assay Type	Sandwich
Detection Range	62.5 pg/ml -4000 pg/ml.
Sensitivity	The minimum detectable dose of human KLK5 is typically less than 15.6 pg/ml. The sensitivity of this assay, or Lower Limit of Detection (LLD) was defined as the lowest protein concentration that could be differentiated from zero. It was determined the mean O.D value of 20 replicates of the zero standard added by their three standard deviations.
Intra-assay Precision	Intra-assay Precision (Precision within an assay): CV%<8%. Three samples of known concentration were tested twenty times on one plate to assess.
Inter-assay Precision	Inter-assay Precision (Precision between assays): CV%<10%. Three samples of known concentration were tested in twenty assays to assess.
Preparation and Storage	Unopened test kits should be stored at 2 to 8 degree C upon receipt. Please refer to pdf manual for further storage instructions.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for KLK5 elisa kit	Principle of the assay: This assay employs the quantitative sandwich enzyme immunoassay technique. Antibody specific for KLK5 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any KLK5 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for KLK5 is added to the wells. After washing, avidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of KLK5 bound in the initial step. The color development is stopped and the intensity of the color is measured.
--	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	117306172
NCBI GeneID	25818
NCBI Accession #	NP_001070959.1
NCBI GenBank Nucleotide #	NM_001077491.1
UniProt Accession #	Q9Y337; Q53ZR3; Q9HBG8
Molecular Weight	32,019 Da
NCBI Official Full Name	kallikrein-5 preproprotein
NCBI Official Synonym Full Names	kallikrein-related peptidase 5
NCBI Official Symbol	KLK5
NCBI Official Synonym Symbols	SCTE; KLKL2; KLK-L2
NCBI Protein Information	kallikrein-5; kallikrein 5; kallikrein-like protein 2; stratum corneum tryptic enzyme
UniProt Protein Name	Kallikrein-5
UniProt Gene Name	KLK5
UniProt Synonym Gene Names	SCTE; UNQ570/PRO1132; KLK-L2
UniProt Entry Name	KLK5_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/klk5-human-elisa-kits/kallikrein-related-peptidase-5/112782>