

Rabbit coagulation factor VII (serum prothrombin conversion accelerator) ELISA Kit | F7 elisa kit

Rabbit Coagulation factor VII, F7 ELISA Kit

CATALOG # AAA112806	REACTIVITY Rabbit
-------------------------------	-----------------------------


PRODUCT SYNONYMS

coagulation factor VII (serum prothrombin conversion accelerator); Rabbit Coagulation factor VII; F7 ELISA Kit; Rabbit Coagulation factor VII (F7) ELISA kit; FVII coagulation protein; coagulation factor VII; eptacog alfa; F7 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Rabbit
Specificity	This assay has high sensitivity and excellent specificity for detection of rabbit F7. No significant cross-reactivity or interference between rabbit F7 and analogues was observed.
Sequence Length	444
Samples	serum, plasma, tissue homogenates.
Assay Type	Competitive (Quantitative)
Detection Range	3.9 ng/ml-250 ng/ml.
Sensitivity	Typically less than 1 ng/ml
Intra-assay Precision	Intra-assay Precision (Precision within an assay): CV%<8%. Three samples of known concentration were tested twenty times on one plate to assess.
Inter-assay Precision	Inter-assay Precision (Precision between assays): CV%<10%. Three samples of known concentration were tested in twenty assays to assess.
Preparation and Storage	Unopened test kits should be stored at 2 to 8 degree C upon receipt. Please refer to pdf manual for further storage instructions.

IMAGES



Standard Curve (Sample)

NCBI AND UNIPROT INFORMATION

NCBI GI #	130495948
NCBI GeneID	100009399
NCBI Accession #	NP_001076148.1

NCBI GenBank Nucleotide #	NM_001082679.1
UniProt Accession #	P98139; P79224
Molecular Weight	49,011 Da
NCBI Official Full Name	coagulation factor VII
NCBI Official Symbol	F7
NCBI Protein Information	coagulation factor VII; serum prothrombin conversion accelerator
UniProt Protein Name	Coagulation factor VII
UniProt Gene Name	F7
UniProt Entry Name	FA7_RABIT

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/f7-rabbit-elisa-kits/coagulation-factor-vii-serum-prothrombin-conversion-accelerator/112806>