

Lysosome-associated membrane glycoprotein 2 (Lamp2) Recombinant Protein | Lamp2 recombinant protein

Recombinant Rat Lysosome-associated membrane glycoprotein 2 (Lamp2)

CATALOG # AAA113507	HOST E Coli or Yeast or Baculovirus or Mammalian Cell
-------------------------------	---

PRODUCT SYNONYMS

Lysosome-associated membrane glycoprotein 2 (Lamp2); Recombinant Rat Lysosome-associated membrane glycoprotein 2 (Lamp2); Recombinant Lysosome-associated membrane glycoprotein 2 (Lamp2); Lysosome-associated membrane glycoprotein 2; LAMP-2; Lysosome-associated membrane protein 2; CD107 antigen-like family member B LGP-110 LGP-96 Lysosomal membrane glycoprotein type B; LGP-B C; Lamp2 recombinant protein

PRODUCT SPECIFICATIONS

Host	E Coli or Yeast or Baculovirus or Mammalian Cell
Purity/Purification	Greater or equal to 85% purity as determined by SDS-PAGE.
Form/Format	Lyophilized or liquid (Format to be determined during the manufacturing process)
Sequence Positions	26-376
Species	Rattus norvegicus (Rat)
Preparation and Storage	Store at -20 degree C. For extended storage, store at -20 or -80 degree C.

NCBI AND UNIPROT INFORMATION

NCBI GI #	126382
NCBI GeneID	24944
UniProt Accession #	P17046
Molecular Weight	45,164 Da
NCBI Official Full Name	Lysosome-associated membrane glycoprotein 2
NCBI Official Synonym Full Names	lysosomal-associated membrane protein 2
NCBI Official Symbol	Lamp2
NCBI Protein Information	lysosome-associated membrane glycoprotein 2; LGP-B; LAMP-2; LGP-96; LGP-110; lysosomal membrane glycoprotein 2; CD107 antigen-like family member B; lysosomal membrane glycoprotein type B; lysosome-associated membrane protein 2
UniProt Protein Name	Lysosome-associated membrane glycoprotein 2
UniProt Gene Name	Lamp2
UniProt Synonym Gene Names	Lamp-2; LAMP-2; Lysosome-associated membrane protein 2; LGP-B
UniProt Entry Name	LAMP2_RAT

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.