

Flagellar hook protein flgE (flgE) Recombinant Protein | flgE recombinant protein

Recombinant Salmonella typhi Flagellar hook protein flgE (flgE)

CATALOG # AAA114336	HOST E Coli or Yeast or Baculovirus or Mammalian Cell
-------------------------------	---


PRODUCT SYNONYMS

Flagellar hook protein flgE (flgE); Recombinant Salmonella typhi Flagellar hook protein flgE (flgE); flgE recombinant protein

PRODUCT SPECIFICATIONS

Host	E Coli or Yeast or Baculovirus or Mammalian Cell
Purity/Purification	Greater or equal to 85% purity as determined by SDS-PAGE.
Form/Format	Lyophilized or liquid (Format to be determined during the manufacturing process)
Sequence Positions	2-403, Full length protein
Sequence	SFSQAVSGLN AAATNLDVIG NNANSATYG FKSGTASFAD MFAGSKVGLG VKVAGITQDF TDGTTTNTGR GLDVAISQNG FFRLVDSNGS VFYSRNGQFK LDENRNLVNM QGMQLTGYP TGPPTIQQG ANPAPITIPN TLMAAKSTTT ASMQINLNST DPVPSKTPFS VSDADSYNKK GTVTYDSQG NAHDMNVYFV KTKDNEWAVY THDSSDPAAT APTTASTTLK FNENGILESG GTVNITTGTI NGATAATFSL SFLNSMQQNT GANNIVATNQ NGYKPGDLVS YQINNDGTVV GNYSNEQEQV LGQIVLANFA NNEGLASQGD NVWAATQASG VALLGTAGSG NFGKLTNGAL EASNVDLSKE LVNMIVAQRN YQSNAQTIKT QDQILNTLVN LR
Sequence Length	402
Species	Salmonella typhi
Preparation and Storage	Store at -20 degrees C. For long-term storage, store at -20 degrees C or -80 degrees C. Store working aliquots at 4 degrees C for up to one week. Repeated freezing and thawing is not recommended.

IMAGES



SDS-PAGE

NCBI AND UNIPROT INFORMATION

NCBI GI #	16760053
NCBI GeneID	1247620
NCBI Accession #	NP_455670.1
NCBI GenBank Nucleotide #	NC_003198.1
UniProt Accession #	P0A1J2; P16322
Molecular Weight	42,210 Da
NCBI Official Full Name	flagellar hook protein FlgE
NCBI Official Symbol	STY1216

NCBI Official Synonym Symbols	flgE
NCBI Protein Information	flagellar hook protein FlgE
UniProt Protein Name	Flagellar hook protein FlgE
UniProt Gene Name	flgE
UniProt Synonym Gene Names	fla FV; flaK

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/flge-recombinant-protein-e-coli-or-yeast-or-baculovirus-or-mammalian-cell/114336>