

ATP-dependent protease ATPase subunit HslU (hslU) Recombinant Protein | hslU recombinant protein

Recombinant Escherichia coli O157:H7 ATP-dependent protease ATPase subunit HslU (hslU)

CATALOG #
AAA114478

HOST
E Coli or Yeast or Baculovirus or Mammalian Cell

PRODUCT SYNONYMS

ATP-dependent protease ATPase subunit HslU (hslU); Recombinant Escherichia coli O157:H7 ATP-dependent protease ATPase subunit HslU (hslU); hslU recombinant protein

PRODUCT SPECIFICATIONS

Host	E Coli or Yeast or Baculovirus or Mammalian Cell
Purity/Purification	Greater or equal to 85% purity as determined by SDS-PAGE.
Form/Format	Lyophilized or liquid (Format to be determined during the manufacturing process)
Sequence Positions	1-443aa; Full Length Protein
Sequence	MSEMPREIV SELDKHIIGQ DNAKRSVAIA LRNRWRRMQL NEELRHEVTP KNILMIGPTG VGKTEIARRL AKLANAPFIK VEATKFTEVG YVGKEVDSII RDLTDAAVKM VRVQAIEKNR YRAEELAEER ILDLVIPPAAK NNWGQTEQQQ EPSAARQAFR KKLREGQLDD KEIEDLAAA PMGVEIMAPP GMEEMTSQLQ SMFQNLGGQK QKARKLKIKD AMKLLIEEEA AKLVNPEELK QDAIDAVEQH GIVFIDEIDK ICKRGESSGP DVSREGVQRD LLPLVEGCTV STKHGMVKTD HILFIASGAF QIAKPSDLIP ELQGRLPPIRV ELQALTTSD F ERILTEPNAS ITVQYKALMA TEGVNIEFTD SGIKRIAEAA WQVNESTENI GARRLHTVLE RLMEEISYDA SDLSGQNITI DADYVSKHLD ALVADEDLSR FIL
Species	Escherichia coli O157:H7
Preparation and Storage	Store at -20 degrees C. For long-term storage, store at -20 degrees C or -80 degrees C. Store working aliquots at 4 degrees C for up to one week. Repeated freezing and thawing is not recommended.

NCBI AND UNIPROT INFORMATION

NCBI GI #	15834112
NCBI GeneID	915031
NCBI Accession #	NP_312885.1
NCBI GenBank Nucleotide #	NC_002695.1
UniProt Accession #	P0A6H6; P32168
Molecular Weight	49,594 Da
NCBI Official Full Name	ATP-dependent protease ATP-binding protein HslU
NCBI Official Symbol	hslU
NCBI Protein Information	ATP-dependent protease ATP-binding protein HslU
UniProt Protein Name	ATP-dependent protease ATPase subunit HslU
UniProt Gene Name	hslU

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/hslu-recombinant-protein-e-coli-or-yeast-or-baculovirus-or-mammalian-cell/114478>