

E3 ubiquitin-protein ligase ICP0 (ICP0) Recombinant Protein | ICP0 recombinant protein

Recombinant Human herpesvirus 1 E3 ubiquitin-protein ligase ICP0 (ICP0), partial

CATALOG #
AAA114921

HOST
E Coli or Yeast or Baculovirus or Mammalian Cell

PRODUCT SYNONYMS

E3 ubiquitin-protein ligase ICP0 (ICP0); Recombinant Human herpesvirus 1 E3 ubiquitin-protein ligase ICP0 (ICP0); partial; ICP0 recombinant protein

PRODUCT SPECIFICATIONS

Host	E Coli or Yeast or Baculovirus or Mammalian Cell
Purity/Purification	Greater or equal to 85% purity as determined by SDS-PAGE.
Form/Format	Lyophilized or liquid (Format to be determined during the manufacturing process)
Sequence Positions	1-300aa; partial
Sequence	MEPRPGASTR RPEGRPQREP APDVWFPCD RDLPDSSDSE AETEVGGRGD ADHHDDDSAS EADSTDTELF ETGLLGPQGV DGGAVSGGSP PREEDPGSCG GAPPREDGGS DEGDVCAVCT DEIAPHLRCD TFPCMRFCI PCMKTWMQLR NTCPLCNAKL VYLIVGVTPS GSFTIPIVN DPQTRMEAAE AVRAGTAVDF IWTGNQRFAP RYLTLGGHTV RALSPHPEP TTDEDDDDLD DADYVPPAPR RTPRAPPRRG AAAPPVTGGA SHAAPQPA RTAPPSAPIG PHGSSNTNTT
Species	Human herpesvirus 1 (strain 17) (HHV-1) (Human herpes simplex virus 1)
Preparation and Storage	Store at -20 degrees C. For long-term storage, store at -20 degrees C or -80 degrees C. Store working aliquots at 4 degrees C for up to one week. Repeated freezing and thawing is not recommended.

NCBI AND UNIPROT INFORMATION

NCBI GI #	820945229
NCBI GeneID	2703389
NCBI Accession #	YP_009137074.1
NCBI GenBank Nucleotide #	NC_001806.2
UniProt Accession #	P08393
Molecular Weight	78,457 Da
NCBI Official Full Name	ubiquitin E3 ligase ICP0
NCBI Official Symbol	RL2
NCBI Protein Information	contains a RING finger; disrupts ND10; proteasome-dependent degradation of several cellular proteins; involved in gene regulation; involved in cellular protein degradation; involved in latency
UniProt Protein Name	E3 ubiquitin-protein ligase ICP0
UniProt Gene Name	ICP0
UniProt Synonym Gene Names	IE110

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

