

Outer membrane protein assembly factor YaeT (yaeT) Recombinant Protein | yaeT recombinant protein

Recombinant *Shigella flexneri* Outer membrane protein assembly factor YaeT (yaeT), partial

CATALOG #
AAA115148

HOST
E Coli or Yeast or Baculovirus or Mammalian Cell

PRODUCT SYNONYMS

Outer membrane protein assembly factor YaeT (yaeT); Recombinant *Shigella flexneri* Outer membrane protein assembly factor YaeT (yaeT); partial; yaeT recombinant protein

PRODUCT SPECIFICATIONS

Host	E Coli or Yeast or Baculovirus or Mammalian Cell
Purity/Purification	Greater or equal to 85% purity as determined by SDS-PAGE.
Form/Format	Lyophilized or liquid (Format to be determined during the manufacturing process)
Sequence Positions	24-421aa; partial
Sequence	FVVKDIHFEG LQRVAVGAAL LSMPVRTGDT VNDEDISNTI RALFATGNFE DVRVLRDGDG LLVQVKERPT IASITFSGNK SVKDDMLKQN LEASGVRVGE SLDRTTIADI EKGLEDFYYS VGKYSASVKA VVTPLPNRV DLKLVFQEGV SAEIQQINIV GNHAFTTDEL ISHFQLRDEV PWWNVVGDRK YQKQKLAGDL ETLRSYYLDR GYARFNIDST QVSLTPDKKG IYVTVNITEG DQYKLSGVEV SGNLAGHSAE IEQLTKIEPG ELYNGTKVTK MEDDIKLLG RYGYAYPRVQ SMPEINDADK TVKLRVNDA GNRFYVRKIR FEGNDTSKDA VLRREMRQME GAWLGSDLVD QGKERLNRLG FFETVDTDTQ RVPGPSDQVD VVYKVKER
Species	<i>Shigella flexneri</i>
Preparation and Storage	Store at -20 degrees C. For long-term storage, store at -20 degrees C or -80 degrees C. Store working aliquots at 4 degrees C for up to one week. Repeated freezing and thawing is not recommended.

NCBI AND UNIPROT INFORMATION

NCBI GI #	24111612
NCBI GeneID	1024444
NCBI Accession #	NP_706122.1
NCBI GenBank Nucleotide #	NC_004337.2
UniProt Accession #	P0A943 ; P39170 ; P39181 ; P77465 ; Q549H3
Molecular Weight	90,553 Da
NCBI Official Full Name	outer membrane protein assembly factor YaeT
NCBI Official Symbol	yaeT
NCBI Protein Information	outer membrane protein assembly factor YaeT
UniProt Protein Name	Outer membrane protein assembly factor BamA
UniProt Gene Name	bamA

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/bama-recombinant-protein-e-coli-or-yeast-or-baculovirus-or-mammalian-cell/115148>