

Nucleoside-specific channel-forming protein tsx (tsx) Recombinant Protein | tsx recombinant protein

Recombinant Escherichia coli Nucleoside-specific channel-forming protein tsx (tsx)

CATALOG # AAA116107	HOST E Coli or Yeast or Baculovirus or Mammalian Cell
-------------------------------	---

PRODUCT SYNONYMS

Nucleoside-specific channel-forming protein tsx (tsx); Recombinant Escherichia coli Nucleoside-specific channel-forming protein tsx (tsx); tsx recombinant protein

PRODUCT SPECIFICATIONS

Host	E Coli or Yeast or Baculovirus or Mammalian Cell
Purity/Purification	Greater or equal to 85% purity as determined by SDS-PAGE.
Form/Format	Lyophilized or liquid (Format to be determined during the manufacturing process)
Sequence Positions	23-294aa; Full Length of Mature Protein
Sequence	AENDKPQYLS DWWHQSVNVV GSYHTRFGPQ IRNDTYLEYE AFAKKDWDFD YGYADAPVFF GGNSDAKGIW NHGSPLFMEI EPRFSIDKLT NTDLSFGPFK EWYFANNYIY DMGRNKDGRQ STWYMGLGTD IDTGLPMSLS MNVYAKYQWQ NYGAANENEW DGYRFKIKYF VPITDLWGGQ LSYIGFTNFD WGS DLGDDSG NAINGIKTRT NNSIASSHIL ALNYDHWYHS VVARYWHDGG QWNDDAELNF GNGNFNVRST GWGGYLVVGY NF
Species	Escherichia coli (strain K12)
Preparation and Storage	Store at -20 degrees C. For long-term storage, store at -20 degrees C or -80 degrees C. Store working aliquots at 4 degrees C for up to one week. Repeated freezing and thawing is not recommended.

NCBI AND UNIPROT INFORMATION

NCBI GI #	16128396
NCBI GeneID	946242
NCBI Accession #	NP_414945.1
NCBI GenBank Nucleotide #	NC_000913.3
UniProt Accession #	P0A927 ; P22786 ; Q2MC15
Molecular Weight	33,589 Da
NCBI Official Full Name	nucleoside channel, receptor of phage T6 and colicin K
NCBI Official Symbol	tsx
NCBI Official Synonym Symbols	ECK0405; JW0401; nupA; T6rec
NCBI Protein Information	nucleoside channel, receptor of phage T6 and colicin K
UniProt Protein Name	Nucleoside-specific channel-forming protein Tsx
UniProt Gene Name	tsx

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

