

Ras-related and estrogen-regulated growth inhibitor-like protein (RERGL) Recombinant Protein | RERGL recombinant protein

Recombinant Human Ras-related and estrogen-regulated growth inhibitor-like protein (RERGL)

CATALOG # AAA116894	HOST E Coli or Yeast or Baculovirus or Mammalian Cell
-------------------------------	---

PRODUCT SYNONYMS

Ras-related and estrogen-regulated growth inhibitor-like protein (RERGL); Recombinant Human Ras-related and estrogen-regulated growth inhibitor-like protein (RERGL); RERGL recombinant protein

PRODUCT SPECIFICATIONS

Host	E Coli or Yeast or Baculovirus or Mammalian Cell
Purity/Purification	Greater or equal to 85% purity as determined by SDS-PAGE.
Form/Format	Lyophilized or liquid (Format to be determined during the manufacturing process)
Sequence Positions	1-205aa; Full Length
Sequence	MSNFLHLKYN EKSVSVTAL TVRFLTKRFI GEYASNFESI YKKHLCLERK QLNLEIYDPC SQTQKAKFSL TSELHWADGF VIVYDISDRS SFAFAKALIY RIREPQTSHC KRAVESAVFL VGNKRDLCVH REVGWEEGQK LALENRCQFC ELSAAEQSLE VEMMFIRIHK DILINFKLKE KRRPSGSKSM AKLINNVFGK RRKSV
Species	Human
Preparation and Storage	Store at -20 degrees C. For long-term storage, store at -20 degrees C or -80 degrees C. Store working aliquots at 4 degrees C for up to one week. Repeated freezing and thawing is not recommended.

NCBI AND UNIPROT INFORMATION

NCBI GI #	555290053
NCBI GeneID	79785
NCBI Accession #	NP_001273130.1
NCBI GenBank Nucleotide #	NM_001286201.1
UniProt Accession #	Q9H628
Molecular Weight	23,855 Da
NCBI Official Full Name	ras-related and estrogen-regulated growth inhibitor-like protein isoform b
NCBI Official Synonym Full Names	RERG like
NCBI Official Symbol	RERGL
NCBI Protein Information	ras-related and estrogen-regulated growth inhibitor-like protein
UniProt Protein Name	Ras-related and estrogen-regulated growth inhibitor-like protein
UniProt Gene Name	RERGL

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

