

Killer cell lectin-like receptor 3 (Klra3) Recombinant Protein | Klra3 recombinant protein

Recombinant Mouse Killer cell lectin-like receptor 3 (Klra3), partial

CATALOG #
AAA116901

HOST
E Coli or Yeast or Baculovirus or Mammalian Cell

PRODUCT SYNONYMS

Killer cell lectin-like receptor 3 (Klra3); Recombinant Mouse Killer cell lectin-like receptor 3 (Klra3); partial; Klra3 recombinant protein

PRODUCT SPECIFICATIONS

Host	E Coli or Yeast or Baculovirus or Mammalian Cell
Purity/Purification	Greater or equal to 85% purity as determined by SDS-PAGE.
Form/Format	Lyophilized or liquid (Format to be determined during the manufacturing process)
Sequence Positions	70-266aa. Partial, provide the complete extracellular domain.
Sequence	QYNQHKQEIN ETLNHHHNCs NMQRAFNLKE EMLTNKSIDC RPSNETLEYI KREQDRWDSK TKTVLDSSRD TGRGVKYWFC YSTKCYFIM NKTTWSGCKA NCQHYSVPIL KIEDEDELKF LQRHVIPENY WIGLSYDKKK KEAWAWIDNGP SKLDMKIRKM NFKSRGCVFL SKARIEDIDC NIPYYCICGK KLDKFPD
Species	Mouse
Preparation and Storage	Store at -20 degrees C. For long-term storage, store at -20 degrees C or -80 degrees C. Store working aliquots at 4 degrees C for up to one week. Repeated freezing and thawing is not recommended.

NCBI AND UNIPROT INFORMATION

NCBI GI #	575771850
NCBI GeneID	16634
NCBI Accession #	NP_001276533.1
NCBI GenBank Nucleotide #	NM_001289604.1
UniProt Accession #	Q64329 ; Q3TWQ1 ; Q3ZB40 ; Q499I4 ; Q61154 ; Q61198 ; Q64257
Molecular Weight	27.6 kDa
NCBI Official Full Name	killer cell lectin-like receptor 3
NCBI Official Synonym Full Names	killer cell lectin-like receptor, subfamily A, member 3
NCBI Official Symbol	Klra3
NCBI Official Synonym Symbols	5E6; Nk2; Nk-2; Ly49c; Nk2.1; NK-2.1
NCBI Protein Information	killer cell lectin-like receptor 3
UniProt Protein Name	Killer cell lectin-like receptor 3
UniProt Gene Name	Klra3
UniProt Synonym Gene Names	Ly-49c; Ly49C; Ly-49c

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

