

Serine/threonine-protein kinase pim-2 (PIM2) Recombinant Protein | PIM2 recombinant protein

Recombinant Human Serine/threonine-protein kinase pim-2 (PIM2)

CATALOG # AAA117105	HOST E Coli or Yeast or Baculovirus or Mammalian Cell
-------------------------------	---

PRODUCT SYNONYMS

Serine/threonine-protein kinase pim-2 (PIM2); Recombinant Human Serine/threonine-protein kinase pim-2 (PIM2); PIM2 recombinant protein

PRODUCT SPECIFICATIONS

Host	E Coli or Yeast or Baculovirus or Mammalian Cell
Purity/Purification	Greater or equal to 85% purity as determined by SDS-PAGE.
Form/Format	Lyophilized or liquid (Format to be determined during the manufacturing process)
Sequence Positions	1-311aa; Full Length
Sequence	MLTKPLQGPP APPGTPTPPP GSKDREAFAE EYRLGPLLGK GGFGTVFAGH RLTDRLQVAI KVIPRNRVLG WSPLSDSVTC PLEVALLWKV GAGGGHPGVI RLLDWFETQE GFMLVLERPL PAQDLFDYIT EKGPLGEGPS RCFFGQVVAA IQHCHSRGVV HRDIKDENIL IDLRGCAKL IDFGSGALLH DEPYTDFDGT RVYSPPEWIS RHQYHALPAT VWSLGILLYD MCVGDIPFER DQEILEAEH FPAHVSPDCC ALIRRCLAPK PSSRPSLEEI LLDPWMQTPA EDVPLNPSKG GPAPLAWSLL P
Species	Human
Preparation and Storage	Store at -20 degrees C. For long-term storage, store at -20 degrees C or -80 degrees C. Store working aliquots at 4 degrees C for up to one week. Repeated freezing and thawing is not recommended.

NCBI AND UNIPROT INFORMATION

NCBI GI #	42821112
NCBI GeneID	11040
NCBI Accession #	NP_006866.2
NCBI GenBank Nucleotide #	NM_006875.3
UniProt Accession #	Q9P1W9 ; Q99739 ; A8K4G6
Molecular Weight	34,190 Da
NCBI Official Full Name	serine/threonine-protein kinase pim-2
NCBI Official Synonym Full Names	Pim-2 proto-oncogene, serine/threonine kinase
NCBI Official Symbol	PIM2
NCBI Protein Information	serine/threonine-protein kinase pim-2
UniProt Protein Name	Serine/threonine-protein kinase pim-2
UniProt Gene Name	PIM2

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

