

Probable vesicular-fusion protein sec17 homolog (B1D1.150, NCU01674) Recombinant Protein | B1D1.150 recombinant protein

Recombinant Neurospora crassa Probable vesicular-fusion protein sec17 homolog (B1D1.150, NCU01674)

CATALOG #
AAA117244

HOST
E Coli or Yeast or Baculovirus or Mammalian Cell

PRODUCT SYNONYMS

Probable vesicular-fusion protein sec17 homolog (B1D1.150; NCU01674); Recombinant Neurospora crassa Probable vesicular-fusion protein sec17 homolog (B1D1.150; B1D1.150 recombinant protein

PRODUCT SPECIFICATIONS

Host	E Coli or Yeast or Baculovirus or Mammalian Cell
Purity/Purification	Greater or equal to 85% purity as determined by SDS-PAGE.
Form/Format	Lyophilized or liquid (Format to be determined during the manufacturing process)
Sequence Positions	1-292aa; Full Length
Sequence	MAVDPRVLQQ EAEKTLASAS KGWGLFGNKE DKYQNAADQY IQAANAFRLQ KSNTEAGKCF EEAAKIFTEK LKEPNDAAANA MLDAFKVYRK DAPDNAVRCV EVAIKQYTMA GNFRRAASHK ENQAEVYENE LQNKPEAIKA YTAAEWYEN DGAVALANKL WLKVADLSAL AGDFFAAIEK FEKVAEASLG NNLMRYSVKE YFLKAGLCSL ATKDMVTAQR NITKYAEKDP SFTGQREYQL LVDLLEASN NNLEMFQDKL AAYDKMSRLD DWKAAVLLQI KNNFEEADNE FS
Species	Neurospora crassa (strain 24698 / 74-OR23-1A / CBS 708.71 / DSM 1257 / FGSC 987)
Preparation and Storage	Store at -20 degrees C. For long-term storage, store at -20 degrees C or -80 degrees C. Store working aliquots at 4 degrees C for up to one week. Repeated freezing and thawing is not recommended.

NCBI AND UNIPROT INFORMATION

NCBI GI #	164428115
NCBI GeneID	3873042
NCBI Accession #	XP_956880.2
NCBI GenBank Nucleotide #	XM_951787.3
UniProt Accession #	Q9P6A5; Q7RY28
Molecular Weight	32,805 Da
NCBI Official Full Name	alpha-soluble NSF attachment protein
NCBI Official Symbol	NCU01674
NCBI Protein Information	alpha-soluble NSF attachment protein
UniProt Protein Name	Probable vesicular-fusion protein sec17 homolog
UniProt Gene Name	B1D1.150

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

