

Isoamyl acetate-hydrolyzing esterase 1 homolog (IAH1) Recombinant Protein | IAH1 recombinant protein

Recombinant Bovine Isoamyl acetate-hydrolyzing esterase 1 homolog (IAH1)

CATALOG # AAA117371	HOST E Coli or Yeast or Baculovirus or Mammalian Cell
-------------------------------	---

PRODUCT SYNONYMS

Isoamyl acetate-hydrolyzing esterase 1 homolog (IAH1); Recombinant Bovine Isoamyl acetate-hydrolyzing esterase 1 homolog (IAH1); IAH1 recombinant protein

PRODUCT SPECIFICATIONS

Host	E Coli or Yeast or Baculovirus or Mammalian Cell
Purity/Purification	Greater or equal to 85% purity as determined by SDS-PAGE.
Form/Format	Lyophilized or liquid (Format to be determined during the manufacturing process)
Sequence Positions	1-249aa; Full Length
Sequence	MALCEAVASG SPLVWPRVLL FGDSITQFSF QQGGWGASLA DMLVRKCDVL NRGFSGYNTR WAKIILPRLV RKGSGLDSPV AVTIFFGAND SALKDENPKQ HVPLEEFVAN LRSMVRYLRS VDVPEGRILIL ITPPPLCEAA WAQECLQQGC KLNRLNSVVG EYARACLQVA QDCGADALDL WSLMQKDGQD FSSYLSDGLH LSPKGNEFVF SHLWPLIEKK VSSLPFLLPY WRDIAEARPE LSLGDDGDH
Species	Bovine
Preparation and Storage	Store at -20 degrees C. For long-term storage, store at -20 degrees C or -80 degrees C. Store working aliquots at 4 degrees C for up to one week. Repeated freezing and thawing is not recommended.

NCBI AND UNIPROT INFORMATION

NCBI GI #	77736461
NCBI GeneID	614320
NCBI Accession #	NP_001029930.1
NCBI GenBank Nucleotide #	NM_001034758.1
UniProt Accession #	Q3SZ16
Molecular Weight	27,579 Da
NCBI Official Full Name	isoamyl acetate-hydrolyzing esterase 1 homolog
NCBI Official Synonym Full Names	isoamyl acetate hydrolyzing esterase 1 (putative)
NCBI Official Symbol	IAH1
NCBI Protein Information	isoamyl acetate-hydrolyzing esterase 1 homolog
UniProt Protein Name	Isoamyl acetate-hydrolyzing esterase 1 homolog
UniProt Gene Name	IAH1

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

