

Alpaca anti-Human Serum Albumin/HSA Nanobody Monoclonal Antibody | anti-HSA antibody

Anti-Human Serum Albumin/HSA Nanobody (ALB8)

CATALOG # AAA120666	HOST Alpaca	REACTIVITY Human	APPLICATIONS ELISA
-------------------------------	-----------------------	----------------------------	------------------------------


PRODUCT SYNONYMS

Serum Albumin/HSA Nanobody; Monoclonal Antibody; Anti-Human Serum Albumin/HSA Nanobody (ALB8); ALB; Albumin; anti-HSA antibody

PRODUCT SPECIFICATIONS

Host	Alpaca
Reactivity	Human
Clonality	Monoclonal
Isotype	VHH-8His-Cys-tag
Clone Number	ALB8
Purity/Purification	>95% by SDS-PAGE Purified by Nickel column.
Form/Format	Liquid; 0.01M PBS, pH7.4.
Concentration	1.5 mg/ml (varies by lot)
Applicable Applications for anti-HSA antibody	ELISA
Expression System	Mammalian cells
Preparation and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4 degree C short term (1-2 weeks). Store at -20 degree C 12 months. Store at -80 degree C long term.

IMAGES



SDS-PAGE

SDS-PAGE for Human Serum Albumin/HSA Nanobody (ALB8)

ADDITIONAL INFORMATION

Product Categories/Family for anti-HSA antibody	Nanoantibodies
---	----------------

NCBI AND UNIPROT INFORMATION

NCBI GI #	4502027
NCBI GeneID	213
NCBI Accession #	NP_000468.1
NCBI GenBank Nucleotide #	NM_000477.5
UniProt Accession #	P02768; O95574; P04277; Q13140; Q645G4; Q68DN5; Q6UXK4; Q86YG0; Q9P157; Q9P1I7; Q9UHS3; Q9UJZ0
NCBI Official Full Name	serum albumin preproprotein
NCBI Official Synonym Full Names	albumin
NCBI Official Symbol	ALB
NCBI Official Synonym Symbols	PRO0883; PRO0903; PRO1341
NCBI Protein Information	serum albumin; albumin (32 AA); albumin (AA 34); growth-inhibiting protein 20; cell growth inhibiting protein 42
UniProt Protein Name	Serum albumin
UniProt Gene Name	ALB
UniProt Entry Name	ALBU_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.