

Human Paraoxonase 2 (PON2) ELISA Kit | PON2 elisa kit

Human Paraoxonase 2 (PON2) ELISA Kit

CATALOG #
AAA120788

REACTIVITY
Human


PRODUCT SYNONYMS

Paraoxonase 2 (PON2); Human Paraoxonase 2 (PON2) ELISA Kit; PON2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	94
Samples	Serum, plasma, cell culture supernates, cell lysates, tissue homogenates
Assay Type	Quantitative Sandwich
Detection Range	0.5ng/ml - 100ng/ml
Sensitivity	0.21ng/ml
Intra-Assay	CV<8%
Inter-Assay	CV<10%
Preparation and Storage	Store at 2 to 8 degree C for 6 months.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PON2 elisa kit

Intended Uses: This sandwich kit is for the accurate quantitative detection of human Paraoxonase 2 (also known as PON2) in serum, plasma, cell culture supernates, cell lysates, tissue homogenates.

Principles of the Assay: This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with human PON2 antibody. PON2 present in the sample is added and binds to antibodies coated on the wells. And then biotinylated human PON2 Antibody is added and binds to PON2 in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated PON2 antibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of human PON2. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm.

NCBI AND UNIPROT INFORMATION

NCBI GI #	34189331
NCBI GeneID	5445
Molecular Weight	39,382 Da
NCBI Official Full Name	PON2 protein, partial
NCBI Official Synonym Full Names	paraoxonase 2
NCBI Official Symbol	PON2
NCBI Protein Information	serum paraoxonase/arylesterase 2
UniProt Protein Name	Serum paraoxonase/arylesterase 2
UniProt Gene Name	PON2
UniProt Synonym Gene Names	PON 2; A-esterase 2

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/paraoxonase-2-pon2/120788>