

Rat Hyperglycosylated Chorionic Gonadotropin (hCG) ELISA Kit | hCG elisa kit

Rat Hyperglycosylated Chorionic Gonadotropin (hCG) ELISA Kit

CATALOG #
AAA120976

REACTIVITY
Rat


PRODUCT SYNONYMS

Hyperglycosylated Chorionic Gonadotropin (hCG); Rat Hyperglycosylated Chorionic Gonadotropin (hCG) ELISA Kit; hCG elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Rat
Assay Type	Quantitative Sandwich
Samples	Serum, plasma, cell culture supernates, Ascites, tissue homogenates or other biological fluids
Detection Range	0.5-150mIU/ml
Sensitivity	0.22mIU/ml
Intra-assay Precision	Intra-Assay Precision (Precision within an assay) Three samples of known concentration were tested on one plate to assess intra-assay precision. Intra-Assay: CV<8%
Inter-assay Precision	Inter-Assay Precision (Precision between assays) Three samples of known concentration were tested in separate assays to assess inter-assay precision. CV(%) = SD/mean x 100. Inter-Assay: CV<10%
Preparation and Storage	Store at 2 to 8 degree C for 6 months.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for hCG elisa kit

Principle of the Assay: This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with Rat hCG antibody. hCG present in the sample is added and binds to antibodies coated on the wells. And then biotinylated Rat hCG Antibody is added and binds to hCG in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated hCG antibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of Rat hCG. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm.

Intended Uses: This sandwich kit is for the accurate quantitative detection of Rat Hyperglycosylated Chorionic Gonadotropin (also known as hCG) in serum, plasma, cell culture supernates, Ascites, tissue homogenates or other biological fluids.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/rat-elisa-kits/hyperglycosylated-chorionic-gonadotropin-hcg/120976>