

Guinea Pig Collagen Type IV, Col IV ELISA Kit | Col IV elisa kit

Guinea Pig Collagen Type IV, Col IV ELISA Kit

CATALOG #
AAA121299

REACTIVITY
Guinea Pig


PRODUCT SYNONYMS

Collagen Type IV; Col IV; Guinea Pig Collagen Type IV; Col IV ELISA Kit; Col IV elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Guinea Pig
Samples	Serum, plasma, cell culture supernates, cell lysates, tissue homogenates
Assay Type	Quantitative Sandwich
Detection Range	2ng/ml - 600ng/ml
Sensitivity	0.97ng/ml
Intra-Assay	CV<8%
Inter-Assay	CV<10%
Preparation and Storage	Store at 2 to 8 degree C for 6 months.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Col IV elisa kit

Intended Uses: This sandwich kit is for the accurate quantitative detection of Guinea Pig Collagen Type IV (also known as Col IV) in serum, plasma, cell culture supernates, cell lysates, tissue homogenates.

Principles of the Assay: This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with Guinea Pig Col IV antibody. Col IV present in the sample is added and binds to antibodies coated on the wells. And then biotinylated Guinea Pig Col IV Antibody is added and binds to Col IV in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated Col IV antibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of Guinea Pig Col IV. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/guinea-pig-pig-elisa-kits/collagen-type-iv-col-iv/121299>