

Porcine adenosine deaminase (ADA) ELISA Kit | ADA elisa kit

Porcine adenosine deaminase (ADA) ELISA Kit

CATALOG #
AAA121383

REACTIVITY
Porcine

PRODUCT SYNONYMS

adenosine deaminase (ADA); Porcine adenosine deaminase (ADA) ELISA Kit; ADA elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Porcine
Sequence Length	363
Samples	Serum, plasma, cell culture supernates, cell lysates, tissue homogenates
Assay Type	Quantitative Sandwich
Detection Range	0.05ng/ml - 20ng/ml
Sensitivity	0.024ng/ml
Intra-Assay	CV<8%
Inter-Assay	CV<10%
Preparation and Storage	Store at 2 to 8 degree C for 6 months.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for ADA elisa kit

Intended Uses: This sandwich kit is for the accurate quantitative detection of Porcine Adenosine Deaminase (also known as ADA) in serum, plasma, cell culture supernates, cell lysates, tissue homogenates.

Principles of the Assay: This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with Porcine ADA antibody. ADA present in the sample is added and binds to antibodies coated on the wells. And then biotinylated Porcine ADA Antibody is added and binds to ADA in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated ADA antibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of Porcine ADA. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm.

NCBI AND UNIPROT INFORMATION

NCBI GI #	28380
NCBI GeneID	100
UniProt Accession #	P00813; Q53F92; Q6LA59
Molecular Weight	40,764 Da
NCBI Official Full Name	adenosine deaminase
NCBI Official Synonym Full Names	adenosine deaminase
NCBI Official Symbol	ADA
NCBI Protein Information	adenosine deaminase
UniProt Protein Name	Adenosine deaminase
UniProt Gene Name	ADA
UniProt Synonym Gene Names	ADA1

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/porcine-elisa-kits/adenosine-deaminase-ada/121383>