

Rabbit 25 Hydroxy Vitamin D3 (HVD3) ELISA Kit | HVD3 elisa kit

Rabbit 25 Hydroxy Vitamin D3 (HVD3) ELISA Kit

CATALOG #
AAA121780

REACTIVITY
Rabbit


PRODUCT SYNONYMS

25 Hydroxy Vitamin D3 (HVD3); Rabbit 25 Hydroxy Vitamin D3 (HVD3) ELISA Kit; HVD3 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Rabbit
Samples	Serum, plasma, cell culture supernates, cell lysates, tissue homogenates
Assay Type	Quantitative Sandwich
Detection Range	0.2nmol/L - 60nmol/L
Sensitivity	0.09nmol/L
Intra-Assay	CV<8%
Inter-Assay	CV<10%
Preparation and Storage	Store at 2 to 8 degree C for 6 months.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for HVD3 elisa kit

Intended Uses: This sandwich kit is for the accurate quantitative detection of Rabbit 25 Hydroxy Vitamin D3 (also known as HVD3) in serum, plasma, cell culture supernates, cell lysates, tissue homogenates.

Principles of the Assay: This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with Rabbit HVD3 antibody. HVD3 present in the sample is added and binds to antibodies coated on the wells. And then biotinylated Rabbit HVD3 Antibody is added and binds to HVD3 in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated HVD3antibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of Rabbit HVD3. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/rabbit-elisa-kits/25-hydroxy-vitamin-d3-hvd3/121780>