

Fish Vitamin D (VD) ELISA Kit | VD elisa kit

Fish Vitamin D (VD) ELISA Kit

CATALOG #
AAA121841

REACTIVITY
Fish


PRODUCT SYNONYMS

Vitamin D (VD); Fish Vitamin D (VD) ELISA Kit; VD elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Fish
Samples	Serum, plasma, cell culture supernates, cell lysates, tissue homogenates
Assay Type	Quantitative Sandwich
Detection Range	0.1?g/ml - 40?g/ml
Sensitivity	0.052?g/ml
Intra-Assay	CV<8%
Inter-Assay	CV<10%
Preparation and Storage	Store at 2 to 8 degree C for 6 months.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for VD elisa kit

Intended Uses: This sandwich kit is for the accurate quantitative detection of Fish Vitamin D (also known as VD) in serum, plasma, cell culture supernates, cell lysates, tissue homogenates.

Principles of the Assay: This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with Fish VD antibody. VD present in the sample is added and binds to antibodies coated on the wells. And then biotinylated Fish VD Antibody is added and binds to VD in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated VD antibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of Fish VD. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/fish-elisa-kits/vitamin-d-vd/121841>