

Fish Lipid Peroxide, LPO ELISA Kit | LPO elisa kit

Fish Lipid Peroxide, LPO ELISA Kit

CATALOG #
AAA121853

REACTIVITY
Fish


PRODUCT SYNONYMS

Lipid Peroxide; LPO; Fish Lipid Peroxide; LPO ELISA Kit; LPO elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Fish
Samples	Serum, plasma, cell culture supernates, cell lysates, tissue homogenates
Assay Type	Quantitative Sandwich
Detection Range	0.5nmol/ml - 150nmol/ml
Sensitivity	0.22nmol/ml
Intra-Assay	CV<8%
Inter-Assay	CV<10%
Preparation and Storage	Store at 2 to 8 degree C for 6 months.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for LPO elisa kit

Intended Uses: This sandwich kit is for the accurate quantitative detection of Fish Lipid Peroxides (also known as LPO) in serum, plasma, cell culture supernates, cell lysates, tissue homogenates.

Principles of the Assay: This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with Fish LPO antibody. LPO present in the sample is added and binds to antibodies coated on the wells. And then biotinylated Fish LPO Antibody is added and binds to LPO in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated LPO antibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of Fish LPO. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/fish-elisa-kits/lipid-peroxide-lpo/121853>