

Fish Hydroxyproline (HYP) ELISA Kit | HYP elisa kit

Fish Hydroxyproline (HYP) ELISA Kit

CATALOG #
AAA121858

REACTIVITY
Fish


PRODUCT SYNONYMS

Hydroxyproline (HYP); Fish Hydroxyproline (HYP) ELISA Kit; HYP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Fish
Samples	Serum, plasma, cell culture supernates, cell lysates, tissue homogenates
Assay Type	Quantitative Sandwich
Detection Range	30ng/ml - 6000ng/ml
Sensitivity	14.33ng/ml
Intra-Assay	CV<8%
Inter-Assay	CV<10%
Preparation and Storage	Store at 2 to 8 degree C for 6 months.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for HYP elisa kit

Intended Uses: This sandwich kit is for the accurate quantitative detection of Fish Hydroxyproline (also known as HYP) in serum, plasma, cell culture supernates, cell lysates, tissue homogenates.

Principles of the Assay: This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with Fish HYP antibody. HYP present in the sample is added and binds to antibodies coated on the wells. And then biotinylated Fish HYP Antibody is added and binds to HYP in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated HYP antibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of Fish HYP. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/fish-elisa-kits/hydroxyproline-hyp/121858>