

Fish Cysteinyl aspartate specific proteinases 3, Caspase-3/CPP32 ELISA Kit

Fish cysteinyl aspartate specific proteinases 3, Caspase-3/CPP32 ELISA Kit

CATALOG #
AAA121890

REACTIVITY
Fish


PRODUCT SYNONYMS

Cysteinyl aspartate specific proteinases 3; Caspase-3/CPP32; Fish cysteinyl aspartate specific proteinases 3; Caspase-3/CPP32 ELISA Kit; Caspase-3/CPP32 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Fish
Samples	Serum, plasma, cell culture supernates, cell lysates, tissue homogenates
Assay Type	Quantitative Sandwich
Detection Range	0.05ng/ml - 10ng/ml
Sensitivity	0.021ng/ml
Intra-Assay	CV<8%
Inter-Assay	CV<10%
Preparation and Storage	Store at 2 to 8 degree C for 6 months.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Cysteinyl aspartate specific proteinases 3, Caspase-3/CPP32 elisa kit

Intended Uses: This sandwich kit is for the accurate quantitative detection of Fish Cysteinyl Aspartate Specific Proteinases 3 (also known as Caspase-3/CPP32) in serum, plasma, cell culture supernates, cell lysates, tissue homogenates.

Principles of the Assay: This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with Fish Caspase-3/CPP32 antibody. Caspase-3/CPP32 present in the sample is added and binds to antibodies coated on the wells. And then biotinylated Fish Caspase-3/CPP32 Antibody is added and binds to Caspase-3/CPP32 in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated Caspase-3/CPP32 antibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of Fish Caspase-3/CPP32. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/fish-elisa-kits/cysteiny-l-aspartate-specific-proteinases-3-caspase-3-cpp32/121890>