

Horse Cardiac troponin I (cTn-I) ELISA Kit | cTn-I elisa kit

Horse Cardiac troponin I (cTn-I) ELISA Kit

CATALOG #
AAA122016

REACTIVITY
Horse


PRODUCT SYNONYMS

Cardiac troponin I (cTn-I); Horse Cardiac troponin I (cTn-I) ELISA Kit; cTn-I elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Horse
Samples	Serum, plasma, cell culture supernates, cell lysates, tissue homogenates
Assay Type	Quantitative Sandwich
Detection Range	5ng/L - 1500ng/L
Sensitivity	2.46ng/L
Intra-Assay	CV<8%
Inter-Assay	CV<10%
Preparation and Storage	Store at 2 to 8 degree C for 6 months.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for cTn-I elisa kit

Intended Uses: This sandwich kit is for the accurate quantitative detection of Horse Cardiac Troponin I (also known as cTn-I) in serum, plasma, cell culture supernates, cell lysates, tissue homogenates.

Principles of the Assay: This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with Horse cTn-I antibody. cTn-I present in the sample is added and binds to antibodies coated on the wells. And then biotinylated Horse cTn-I Antibody is added and binds to cTn-I in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated cTn-I antibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of Horse cTn-I. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/horse-elisa-kits/cardiac-troponin-i-ctn-i/122016>