

Horse major acute phase protein (MAP) ELISA Kit | MAP elisa kit

Horse major acute phase protein (MAP) ELISA Kit

CATALOG #
AAA122018

REACTIVITY
Horse


PRODUCT SYNONYMS

major acute phase protein (MAP); Horse major acute phase protein (MAP) ELISA Kit; MAP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Horse
Samples	Serum, plasma, cell culture supernates, cell lysates, tissue homogenates
Assay Type	Quantitative Sandwich
Detection Range	0.05?g/ml - 10?g/ml`
Sensitivity	0.021?g/ml
Intra-Assay	CV<8%
Inter-Assay	CV<10%
Preparation and Storage	Store at 2 to 8 degree C for 6 months.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for MAP elisa kit

Intended Uses: This sandwich kit is for the accurate quantitative detection of Horse Major Acute Phase Protein (also known as MAP) in serum, plasma, cell culture supernates, cell lysates, tissue homogenates.

Principles of the Assay: This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with Horse MAP antibody. MAP present in the sample is added and binds to antibodies coated on the wells. And then biotinylated Horse MAP Antibody is added and binds to MAP in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated MAP antibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of Horse MAP. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/horse-elisa-kits/major-acute-phase-protein-map/122018>