

Horse Dermatan Sulfate (DS) ELISA Kit | DS elisa kit

Horse Dermatan Sulfate (DS) ELISA Kit

CATALOG #
AAA122022

REACTIVITY
Horse


PRODUCT SYNONYMS

Dermatan Sulfate (DS); Horse Dermatan Sulfate (DS) ELISA Kit; DS elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Horse
Samples	Serum, plasma, cell culture supernates, cell lysates, tissue homogenates
Assay Type	Quantitative Sandwich
Detection Range	0.05ug/ml - 15ug/ml
Sensitivity	0.023ug/ml
Intra-Assay	CV<8%
Inter-Assay	CV<10%
Preparation and Storage	Store at 2 to 8 degree C for 6 months.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for DS elisa kit

Intended Uses: This sandwich kit is for the accurate quantitative detection of Horse Dermatan Sulfate (also known as DS) in serum, plasma, cell culture supernates, cell lysates, tissue homogenates.

Principles of the Assay: This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with Horse DS antibody. DS present in the sample is added and binds to antibodies coated on the wells. And then biotinylated Horse DS Antibody is added and binds to DS in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated DS antibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of Horse DS. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/horse-elisa-kits/dermatan-sulfate-ds/122022>