

Bovine Bone Alkaline Phosphatase (BALP) ELISA Kit | BALP elisa kit

Bovine Bone Alkaline Phosphatase (BALP) ELISA Kit

CATALOG #
AAA122186

REACTIVITY
Bovine


PRODUCT SYNONYMS

Bone Alkaline Phosphatase (BALP); Bovine Bone Alkaline Phosphatase (BALP) ELISA Kit; BALP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Bovine
Samples	Serum, plasma, cell culture supernates, cell lysates, tissue homogenates
Assay Type	Quantitative Sandwich
Detection Range	0.5ng/ml - 200ng/ml
Sensitivity	0.27ng/ml
Intra-Assay	CV<8%
Inter-Assay	CV<10%
Preparation and Storage	Store at 2 to 8 degree C for 6 months.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for BALP elisa kit

Intended Uses: This sandwich kit is for the accurate quantitative detection of Bovine Bone alkaline phosphatase (also known as BALP) in serum, plasma, cell culture supernates, cell lysates, tissue homogenates.

Principles of the Assay: This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with Bovine BALP antibody. BALP present in the sample is added and binds to antibodies coated on the wells. And then biotinylated Bovine BALP Antibody is added and binds to BALP in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated BALPantibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of Bovine BALP. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/bovine-ovine-elisa-kits/bone-alkaline-phosphatase-balp/122186>