

Bovine Fetuin B, FETU-B ELISA Kit | FETU-B elisa kit

Bovine Fetuin B, FETU-B ELISA Kit

CATALOG #
AAA122292

REACTIVITY
Bovine


PRODUCT SYNONYMS

Fetuin B; FETU-B; Bovine Fetuin B; FETU-B ELISA Kit; FETU-B elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Bovine
Samples	Serum, plasma, cell culture supernates, cell lysates, tissue homogenates
Assay Type	Quantitative Sandwich
Detection Range	5ng/ml - 2000ng/ml
Sensitivity	2.48ng/ml
Intra-Assay	CV<8%
Inter-Assay	CV<10%
Preparation and Storage	Store at 2 to 8 degree C for 6 months.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for FETU-B elisa kit

Intended Uses: This sandwich kit is for the accurate quantitative detection of Bovine Fetuin B (also known as FETU-B) in serum, plasma, cell culture supernates, cell lysates, tissue homogenates.

Principles of the Assay: This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with Bovine FETU-B antibody. FETU-B present in the sample is added and binds to antibodies coated on the wells. And then biotinylated Bovine FETU-B Antibody is added and binds to FETU-B in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated FETU-B antibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of Bovine FETU-B. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/bovine-ovine-elisa-kits/fetuin-b-fetu-b/122292>