

Canine Resistin (Retn) ELISA Kit | Retn elisa kit

Canine Resistin (Retn) ELISA Kit

CATALOG #
AAA122371

REACTIVITY
Canine


PRODUCT SYNONYMS

Resistin (Retn); Canine Resistin (Retn) ELISA Kit; Retn elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Sequence Length	108
Samples	Serum, plasma, cell culture supernates, cell lysates, tissue homogenates
Assay Type	Quantitative Sandwich
Detection Range	5ng/L - 1500ng/L
Sensitivity	2.44ng/L
Intra-Assay	CV<8%
Inter-Assay	CV<10%
Preparation and Storage	Store at 2 to 8 degree C for 6 months.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Retn elisa kit

Intended Uses: This sandwich kit is for the accurate quantitative detection of Canine Resistin (also known as Retn) in serum, plasma, cell culture supernates, cell lysates, tissue homogenates.

Principles of the Assay: This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with Canine Retn antibody. Retn present in the sample is added and binds to antibodies coated on the wells. And then biotinylated Canine Retn Antibody is added and binds to Retn in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated Retn antibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of Canine Retn. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm.

NCBI AND UNIPROT INFORMATION

NCBI GI #	301129180
NCBI GeneID	56729
NCBI Accession #	NP_001180303.1
NCBI GenBank Nucleotide #	NM_001193374.1
UniProt Accession #	Q9HD89; Q540D9; Q76B53; D6W649
Molecular Weight	8,686 Da
NCBI Official Full Name	resistin
NCBI Official Synonym Full Names	resistin
NCBI Official Symbol	RETN
NCBI Official Synonym Symbols	ADSF; RSTN; XCP1; FIZZ3; RETN1
NCBI Protein Information	resistin
UniProt Protein Name	Resistin
UniProt Gene Name	RETN
UniProt Synonym Gene Names	FIZZ3; HXCP1; RSTN; ADSF

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/canine-elisa-kits/resistin-retn/122371>